

### 1. (A)

### 1. (B)

## English

#### 1. LOOSEN STRAPS

**(A)** Release the two over-center buckles by pinching the two gray buttons on each of the buckles.

**(B)** Loosen the two remaining straps.

#### 2. POSITION THE BRACE

With the patient in a seated position and leg bent at 90°, place the brace on the leg aligning the center of the hinge with the top of the knee cap. Note: Make sure to push the brace posterior so that the hinges run along the mid-line of the leg.

#### 3. SECURE THE TOP CALF STRAP

Wrap the top calf strap around the calf as high as possible and loop through the D-ring. Pull the tab until a snug and comfortable fit is achieved and secure.

#### 4. SECURE THE LOWER THIGH STRAP

Wrap the lower thigh strap around the thigh and loop through the D-ring. Pull tab until a snug and comfortable fit is achieved and secure.

#### 5. ADJUST THE CRESCENT FINGER PAD

Remove the Crescent Finger Pad from underneath the Y-Pad. Relocate the pad to the desired position along the edge of the patient’s knee cap.

#### 6. SECURE LOWER PATELLO-FORCE STRAP

**(A)** With the patient’s knee at 20° of flexion, wrap the lower Patello-Force Strap around the leg and secure the over-center buckle.

**(B)** Adjust the strap using the fastener at the strap end until a snug and comfortable fit is achieved.

#### 7. SECURE UPPER PATELLO-FORCE STRAP


Repeat Step 6 for the upper Patello-Force Strap.

#### 8. ADJUST THE PATELLO-FORCE STRAP

Use the fasteners for additional adjustment.

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