

Bunion Boot

Instructions for Use CP020191 Rev H 01/24


- 1. Apply the Calf Cuff**
- Open the hook-to-pile closure on the calf cuff.
 - Wrap the cuff around the leg making sure the top of the cuff is in line with the arms of the boot.
 - If necessary, trim the cuff.

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
- 2. Make an Auto-mold Impression**
- Place the foot between the uprights. The ankle bones should be aligned with the uprights and foot centered between the uprights.
 - Stand for 30 seconds to make an impression of the foot in the auto-mold.

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
- 3. Position the Great Toe**
- Trim the toe strap from the lateral side.
 - Reposition the toe pad.
 - Wrap the strap around the toe and secure medially. Trim excess strap if necessary.

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- 4. Secure the Forefoot**
- Fold the forefoot pad over the top of the foot.
 - Starting at the toes, feed the strap through the D-ring. Pull to the desired tension and secure.
 - If necessary, trim excess strap.
 - Repeat for second forefoot strap.

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

**BUNION SURGERY WITHOUT CRUTCHES.
LESS TIME, LESS PAIN, MORE CONVENIENCE.**

5. Secure the Uprights

- Align the uprights with the sides of the leg parallel to the leg bones. When the uprights are aligned, remove the plastic from the uprights.
- Press the uprights against the side of the leg to engage the hook material to the cuff.
- Save the plastic covers for reapplication after washing the cuff.

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
6. Secure the Leg Straps

- Beginning at the ankle, wrap both strap ends firmly forward to engage the straps with the hook material on the outside of the uprights.
- Loop the end of the strap through the D-ring and pull back to tension. Press hook end to strap to engage.
- Repeat process for the remaining leg strap.

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
7

7. Inflating the Air Bladder

After applying the boot using the application instructions provided, locate the hand air pump provided with the boot. Place the clear end of pump into the inflation nozzle located on the right side of boot. Press the clear end of pump into the nozzle firmly to prevent any of the air from escaping. Now squeeze the pump firmly two (2) full times and check the patient's ankle to see how the air bladder is fitting.

The process of inflating the bladder may be done by the patient so that the patient controls the pressures in the bladder. Most patients will only use 2-3 full squeezes of air from the pump. **WARNING:** Do not over inflate the air bladder! Over inflating the bladder may cause the bladder to burst and may stop blood circulation to the patient's foot. Air bladders are designed to take no more than 5 full squeezes of air. Over inflating the air bladder will void its warranty.

Deflating the Air Bladder

To deflate the bladder while in the boot, place the black end of pump into the nozzle and squeeze until the air bladder is flat and all the air has been removed.

NOTE: The air bladder and pad can be adjusted to fit a range of ankle and heel sizes. Simply slide your fingers between the air bladder and the hook on the inside of the boot to peel it away. Now slide the air bladder forward for smaller ankles or slide it back for bigger ankles. Slight tension on the back of the heel is needed before inflating the air bladder.

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This device is offered for sale by or on the order of a physician or other qualified medical professional.

Indications: The Bledsoe Bunion Boot is indicated for post-operative Hallux Abducto Valgus (HAV) surgery when a strict joint position must be maintained. It is also indicated for the treatment of foot indications that would benefit from a more even distribution of body weight over the sole of the foot.

Contraindications: The Bledsoe Bunion Boot is contraindicated for unstable fractures and when the foot width will not fit within one size range larger than a normal size shell.

Caution: Activity modification and proper rehabilitation are an essential part of a safe treatment program. If any additional pain or other symptoms occur while using this device, discontinue use and seek medical attention.

Caution: This device is not intended to prevent injury, or reduce or eliminate the risk of re-injury.

Caution: The white foam layer on this boot may exhibit brittle characteristics at temperatures below 50° F (10° C). Please do not handle the foam while it is cold. Make sure the entire boot is allowed to warm to room temperature before use. This device has been developed for the indications listed above, and may be damaged by unusual activities such as sports, falls, or other accidents. Do not leave this device in a car on a hot day, or in any other place where the temperature may exceed 150°F (65°C) for any extended period of time. To avoid the risk of accident and possible injury to the patient, discontinue use if such an incident has occurred.

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© 2024 Breg, Inc. All Rights Reserved
Carlsbad, CA 92008 U.S.A.