English

5. SET RANGE-OF-MOTION

Depress extension ROM button and slide to desired setting. Depress flexion ROM button and slide to desired setting. Adjust opposite to the same settings of first hinge. NOTE: Settings automatically lock in place. Refer to standard hinge settings instruction sheet enclosed. Caution: The range-of-motion stops on both the medial and lateral hinges must be set the same.

6. LOCK HINGE

Lock hinge at a desired -10°, 0°, 10°, 20°, 30°, and 40° degrees. Bend leg to desired setting and engage lock on each hinge. This will lock the brace in place independent of ROM settings. Make sure ROM settings do not interfere with desired lock setting.

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INDICATIONS: Indicated for immediate application following injury or operative procedures to reduce pain and swelling and aid in healing.

CAUTION: Activity modification and proper rehabilitation are an essential part of a safe treatment program. If any additional pain or other symptoms occur while using this device, discontinue use and seek medical attention.

CAUTION: This device is not intended to prevent injury, or reduce or eliminate the risk of re-injury.

CAUTION: The range-of-motion stops on both the

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medial and lateral hinges must be set the same.

CAUTION: This device has been developed for the indications listed above, and may be damaged by unusual activities such as sports, falls, or other accidents, or exposure of the hinges to debris or liquids. To avoid the risk of accident and possible injury to the patient, discontinue use if such an incident has occurred.

Caution: Do not leave this device in a car on a hot day, or in any other place where the temperature may exceed 150F (65 C.) for any extended period of time.

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
English

1. RELEASE ALL FOUR BUCKLES

Pinch the two buttons on the buckle to release it. Place the brace beside the leg aligning the hinge with the joint line.

2. SET TIBIAL AND THIGH LENGTH

Slide the locking tab down to free the adjustment button. Slide extension arms to desired length and slide locking tab back in place to lock extension arm.

3. APPLY BRACE

(A) G3 FULL FOAM: Remove hinge arm that doesn’t have buckle bases and let it hang free from other hinge arm side. With the leg extended slide the brace underneath the leg with closure at the front. Open both cuffs, wrap cuffs around leg, trim cuffs to desired length, and close cuffs with hook attachment. Make sure the center of the hinge is aligned with the center of the knee. Align the arms of the brace along the midline of the leg. The hinges should be even and parallel with the knee joint. Attach 1” gastroc strap. Wrap around and pull through D-ring then attach to double sided hook and trim off excess strap. Then reattach opposite hinge with center of hinge over center of knee joint.

(B) G3: With the leg extended, slide the brace underneath the leg with closure at the front. Make sure the center of the hinge is aligned with the center of the knee. Align the arms of the brace along the midline of the leg. The hinges should be even and parallel with the knee joint. Attach 1” gastroc strap. Wrap around and pull through D-ring then attach to double sided hook and trim off excess strap.

4. SECURE STRAPS

On arm opposite the over-center buckle, pull anterior strap away from leg disengaging hook allowing strap movement. Reattach over center buckles and remove small piece of hook from strap end. Pull strap to desired tension. Reattach hook end to secure strap and trim off excess length. Repeat process for remaining straps.

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
Continued on reverse side

NOTICE TO USERS
If any serious incident has occurred in relation to the device, report to Breg and the competent authority of the Member State in which the user and/or patient is established.

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© 2023 Breg, Inc. All rights reserved.