**ROADRUNNER**  
da infilare

**SHORTRUNNER**  
da infilare

**ROADRUNNER**  
Wrap-On  
de Envolver  
Enveloppement  
Wickel-Modell  
a fascia

**SHORTRUNNER**  
Wrap-On  
de Envolver  
Enveloppement  
Wickel-Modell  
a fascia

## Italiano

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**CAUTION:** UNITED STATES FEDERAL LAW RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A LICENSED HEALTH CARE PRACTITIONER.

**CAUTION:** FOR SINGLE PATIENT USE ONLY.

**CAUTION:** IF YOU EXPERIENCE INCREASED PAIN, SWELLING, OR ANY ADVERSE REACTIONS WHILE USING THIS PRODUCT, IMMEDIATELY CONSULT YOUR MEDICAL PROFESSIONAL.

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



**WARNING:** CAREFULLY READ USE INSTRUCTIONS AND WARNINGS PRIOR TO OPERATION.

**WARNING:** THIS DEVICE WILL NOT PREVENT OR REDUCE ALL INJURIES. PROPER REHABILITATION AND ACTIVITY MODIFICATION ARE ALSO AN ESSENTIAL PART OF A SAFE TREATMENT PROGRAM.

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


©2023 Breg, Inc. All rights reserved.  
Roadrunner and Shortrunner are trademarks and AirMesh is a registered trademark of Breg, Inc.

EC REP  
E/U authorized representative  
MDSS GmbH  
Schiffgraben 41  
D-30175 Hannover  
Germany

AW-1.02070 Rev N 07/23

## English

### APPLICATION INSTRUCTIONS FOR KNEE BRACES:

1. The wrap-on brace is applied with the closures in the front of the leg and the kneecap cutout centered over the knee. Secure the closures firmly.
2. The slip-on brace is applied by pulling the brace on with the kneecap cutout centered over the knee.
3. Both the wrap-on and pull-on braces have tightening straps. The straps should be tightened firmly, starting with the two straps closest to the knee.

### ADJUSTMENT INSTRUCTIONS Polycentric hinge with stops

4. Both flexion and extension stops are supplied with the brace. The brace comes with a 10° extension stop already installed. The stops are labeled in degrees as follows;

Extension Stops: 0°, 10°, 20°, 30° & 40°  
Flexion Stops: 45°, 60°, 75° & 90° (no flexion stops are installed upon shipment)

To change the stops, first, completely remove the screw as shown. Then remove the stop. Insert the new stop, align with the holes in the hinge, then reinsert the screw. Using a Phillips head screw driver, ensure that the screws are seated in the recess and completely tightened. Both medial (inside) and lateral (outside) hinge stops must be identical in degree.

*NOTE: Due to soft tissue compression, we recommend the use of extension stops at least 10° greater than the desired limit. (e.g., if you desire a 10° extension limit, use the 20° extension stop).*

### TIPS ON USE AND CARE OF YOUR KNEE BRACE:

- After exposure to salt water or dirt, thoroughly rinse the brace with fresh water and allow to air dry.
- A dry lubricant, such as teflon spray, may be used on the hinge occasionally. Do not use wet lubricants.
- The hinge assemblies may be removed to allow easier washing of the sleeve and straps.
- The side bars may be bent to accommodate individual legs. Bend in equal amounts above and below the hinge.
- Periodically inspect the sleeve, straps and hinges for excessive wear. Be sure all screws are tight and stops in place.

### NOTICE TO USERS

If any serious incident has occurred in relation to the device, report to Breg and the competent authority of the Member State in which the user and/or patient is established.

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