X2K®

FITTING INSTRUCTIONS

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

B

⚠ WARNINGS

WARNING: CAREFULLY READ FITTING INSTRUCTIONS AND WARNINGS PRIOR TO USE.

WARNING: IF YOU EXPERIENCE INCREASED PAIN, SWELLING, SKIN IRRITATION, OR ANY ADVERSE REACTIONS WHILE USING THIS PRODUCT, IMMEDIATELY CONSULT YOUR MEDICAL PROFESSIONAL.

WARNING: THIS DEVICE WILL NOT PREVENT OR REDUCE ALL INJURIES. PROPER REHABILITATION AND ACTIVITY MODIFICATION ARE ALSO ESSENTIAL PARTS OF A SAFE TREATMENT PROGRAM. CONSULT YOUR LICENSED HEALTHCARE PROFESSIONAL SHOULD YOU HAVE ANY QUESTIONS REGARDING THE SAFE APPLICATION OF THIS PRODUCT RELATED TO YOUR MEDICAL HISTORY.

CAUTION: CARE AND CLEANING OF THIS PRODUCT IS ESSENTIAL TO ITS CONTINUING STRENGTH AND PERFORMANCE. FOLLOW REMOVAL AND CLEANING PROCESS DESCRIBED BELOW.

CAUTION: WHEN USING BRACES BILATERALLY, MODIFY ACTIVITY AS NECESSARY TO KEEP MEDIAL HINGES FROM CONTACTING EACH OTHER.

CAUTION: FOR SINGLE PATIENT USE ONLY.

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A APPLY BRACE TO LEG

- 1 Unfasten all 5 straps and fold hook/loop tabs back against straps. **Be sure to unfasten front strap above kneecap.**
- 2 Apply brace while sitting with knee bent approximately 30°.
- 3 Position brace so that middle of hinges is even with top of kneecap. This positions brace high, it will settle into correct position with activity.

B FASTEN STRAPS IN ORDER BY FOLLOWING NUMBERS ON STRAP ENDS

- 1 Fasten strap just below knee.
- 2 Fasten rear portion of strap just above knee, ensuring that calf cuff lifts off of shin at least 1" (see illustration above).
- 3 Fasten top strap.
- 4 Fasten bottom strap, pulling cuff to shin.
- 5 Stand up and fasten front portion of strap just above knee.

Tighten straps snugly, but comfortably. After 15 minutes of activity, loosen, then re-tighten all straps in proper order. This is required because of tissue swelling and compression.

CARE AND CLEANING

Straps and pads must be removed and cleaned regularly depending upon your activity level and perspiration (see illustrations D through G.) Hand wash straps and pads in cold water with mild detergent, rinse, AIR DRY. Occasionally apply a dry lubricant to hinge especially after use in water.

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
C HINGE RANGE OF MOTION ADJUSTMENT

Brace is shipped with 10° hinge extension stops installed with additional flexion and extension stops in box. Stops are labeled in degrees:

Extension stops: 0°, 10°, 20°, 30° & 40°
Flexion stops: 45°, 60°, 75° & 90°
(brace is shipped with no flexion stops installed)

To change hinge stops, use screwdriver to remove screws. Remove existing stop, insert new stop, replace screw. Be sure that screws are seated and snug. Both medial (inside) hinge and lateral (outside) hinge stops must be identical in degrees.

D CONDYLE PADS

Brace ships with thick and thin condyle pads which can be interchanged to enhance comfort and fit.

E TOP AND BOTTOM STRAP REMOVAL

Slide strap tabs forward to remove. Loosen straps if necessary.

F REMOVAL OF STRAPS JUST ABOVE AND BELOW KNEE


Locate both strap mounting clips on straps. Unlock clips by depressing red button with ballpoint pen, then slide clips toward front of brace. Remove from post.

G THIGH AND CALF PAD REMOVAL

Each pad is held in place by two snaps, pull to remove.

NOTICE TO USERS

If any serious incident has occurred in relation to the device, report to Breg and the competent authority of the Member State in which the user and/or patient is established.

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© 2023 Breg, Inc.