# X2K® PTO® FITTING INSTRUCTIONS

ANLEITUNG ZUM ANLEGEN  
ISTRUZIONI PER L'USO  
MISE EN PLACE  
INSTRUCCIONES DE AJUSTE



Breg, Inc.  
2382 Faraday Avenue, Suite 300  
Carlsbad, CA 92008 U.S.A.  
P: 800-321-0607  
F: 800-329-2734  
www.breg.com  
© 2023 Breg, Inc.  
AW-1.15360 Rev H 10/23

EC REP

E/U authorized representative  
MDSS GmbH  
Schiffgraben 41  
D-30175 Hannover  
Germany

## ⚠️ WARNINGS

**WARNING:** CAREFULLY READ FITTING INSTRUCTIONS AND WARNINGS PRIOR TO USE.

**WARNING:** IF YOU EXPERIENCE INCREASED PAIN, SWELLING, SKIN IRRITATION, OR ANY ADVERSE REACTIONS WHILE USING THIS PRODUCT, IMMEDIATELY CONSULT YOUR MEDICAL PROFESSIONAL.

**WARNING:** THIS DEVICE WILL NOT PREVENT OR REDUCE ALL INJURIES. PROPER REHABILITATION AND ACTIVITY MODIFICATION ARE ALSO ESSENTIAL PARTS OF A SAFE TREATMENT PROGRAM. CONSULT YOUR LICENSED HEALTHCARE PROFESSIONAL SHOULD YOU HAVE ANY QUESTIONS REGARDING THE SAFE APPLICATION OF THIS PRODUCT RELATED TO YOUR MEDICAL HISTORY.

**CAUTION:** CARE AND CLEANING OF THIS PRODUCT IS ESSENTIAL TO ITS CONTINUING STRENGTH AND PERFORMANCE. FOLLOW REMOVAL AND CLEANING PROCESS DESCRIBED BELOW.

**CAUTION:** WHEN USING BRACES BILATERALLY, MODIFY ACTIVITY AS NECESSARY TO KEEP MEDIAL HINGES FROM CONTACTING EACH OTHER.

**CAUTION:** FOR SINGLE PATIENT USE ONLY.

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## A APPLY BRACE TO LEG

- 1 Unfasten all straps including Patella straps, and fold hook/loop tabs back against straps.
- 2 Apply brace while sitting with knee bent approximately 30°.
- 3 Position brace so that middle of hinges is even with top of kneecap. This positions brace high; it will settle into correct position with activity.

## B FASTEN STRAPS IN ORDER BY FOLLOWING NUMBERS ON STRAP ENDS

- 1 Fasten strap just below knee.
- 2 Fasten strap just above knee, ensuring that calf cuff lifts off of shin at least 1" (see illustration above).
- 3 Fasten top strap.
- 4 Fasten bottom strap, pulling cuff to shin.

## 5 PATELLA PLATE INITIAL FITTING

The Patella Plate position should be adjusted so that the plate is adjacent to the Patella. Remove the Patella Plate from both straps and position. Reattach both straps.

## 6 ADJUST THE PATELLA PLATE STRAPS

Bend the knee approximately 70° and secure the Patella Plate by attaching both top and bottom Patella straps securely.

Strap length can be adjusted. Remove the Tab at the end of the strap. Using a pair of scissors, trim the strap to length then replace the Tab.

Tighten the straps snugly, but comfortably. After 15 minutes of activity, loosen, then retighten all straps in proper order. This is required because of tissue swelling and compression.

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
### C ADJUST THE LATERAL CONDYLE PAD

With brace on, walk 5 minutes to position brace properly and ensure that lateral condyle pad is applying lateral force to knee. If more lateral force is required to relieve symptoms turn upper and lower adjustment screws on side of hinge clockwise. This will move condyle pad inward (toward knee), applying more lateral force to knee.

### D ADJUST TOP AND BOTTOM OF CONDYLE PAD

Upper and lower adjustment screws can be adjusted differently to accommodate various knee contours.

For more comfort, different adjustments of upper and lower screws will distribute lateral force more evenly.

**Important Note:** When wearing brace for an extended period you may experience discomfort caused by pressure from lateral condyle pad. Back off adjustment screws to relieve pressure. Remove brace during periods of inactivity.

### E HINGE RANGE OF MOTION ADJUSTMENT

Brace is shipped with 10° hinge extension stops installed with additional flexion and extension stops in box. Stops are labeled in degrees:

Extension stops: 0°, 10°, 20°, 30° & 40°  
Flexion stops: 45°, 60°, 75° & 90°  
(brace is shipped with no flexion stops installed)

To change hinge stops, use screwdriver to remove screws. Remove existing stop, insert new stop, replace screw. Be sure that screws are seated and snug. Both medial (inside) hinge and lateral (outside) hinge stops must be identical in degrees.

### F TOP AND BOTTOM STRAP REMOVAL

Slide strap tabs forward to remove. Loosen straps if necessary.

### G REMOVAL OF STRAPS JUST ABOVE AND BELOW KNEE

Locate both strap mounting clips on straps. Unlock clips by depressing red button with ballpoint pen, then slide clips toward front of brace. Remove from post.

### THIGH AND CALF PAD REMOVAL


Each pad is held in place by two snaps, pull to remove.

### CARE AND CLEANING

Straps and pads must be removed and cleaned regularly depending upon your activity level and perspiration (see illustrations F and G.) Hand wash straps and pads in cold water with mild detergent, rinse, AIR DRY. Occasionally apply a dry lubricant to hinge especially after use in water.

### NOTICE TO USERS

If any serious incident has occurred in relation to the device, report to Breg and the competent authority of the Member State in which the user and/or patient is established.

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