

Z-13 Knee Brace

CP020252 Rev L 11/23

ENGLISH

1. Unfasten Buckles

Release all four straps on the back side of the brace by squeezing the gray levers on either side of the buckle and withdrawing the over-center buckle. Loosen the straps. (A)

2. Check the Fit

With the knee bent at 45°, lay the brace on the leg and align the hinge. The thigh and calf shell shapes should fit the front of the leg. (B)

If the brace does not fit properly, do not proceed; the brace may need to be adjusted or the leg re-measured.

3. Position the Brace

With the brace still in the prefit position, locate the top of the knee cap. Position the center of the hinge on an imaginary line that goes across the top of the knee cap. Push the hinges posteriorly until the hinges are in the center of the knee joint when viewed from the side.

4. Anchor the Brace

Wrap the upper rear calf strap around the calf as high as possible over the calf muscle. (C)

All straps have additional length adjustment at the opposite end of each strap. See step 7.

5. Fasten Rear Thigh Straps

Following the same instructions as in Step 4, fasten the lower and then upper rear thigh strap.

6. Adjust Front Tibia Strap

Tighten the front tibia strap. (D) If the strap overlaps too far, see step 7 for adjustment instructions.

Center the strap pad.

7. Adjust Strap Length

All straps have additional length adjustment at the end opposite the tab. Release the strap and adjust the length of the strap using the hook and pile fastener. (E)

Special Adjustment for Anterior Cruciate Ligament Injuries


In order to provide increased resistance against tibial subluxation (forward movement of the lower leg relative to the upper leg) when the leg is extended or to reduce force on a reconstructed ACL, it may be desirable to pretension the straps to provide a "shear force" across the knee by increasing the tension on the front calf strap on the tibia and the lower rear thigh strap on the back of the thigh. All straps for ACL patients should be adjusted with the leg bent at 60°-80°. Loosen the front thigh strap and increase the tension on the lower rear thigh strap and, if necessary, the front shin strap to provide a pre-loaded strap force that helps prevent anterior tibial subluxation.

Special Adjustment for Posterior Cruciate Ligament Injuries

In order to provide increased resistance against posterior tibial subluxation (backward movement of the lower leg relative to the upper leg) as the leg bends or to reduce force on a PCL reconstruction, it may be desirable to pretension the straps to provide a "shear force" across the knee by increasing the tension on the upper rear calf strap and the front thigh strap. All straps for PCL patients should be adjusted with the leg in full extension. Loosen the lower rear thigh strap and the front calf strap and increase the tension on the front thigh strap and the upper rear calf strap to provide a pre-loaded shear force that prevents posterior tibial subluxation.

NOTE: If needed, apply the included spacers behind the condyle pads to achieve proper condular pressure on the medial and/or lateral side of the knee.

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