### Hinge range of motion adjustment

**A.** Brace is shipped with 10° hinge extension stops installed with additional flexion and extension stops in box. Stops are labeled in degrees:

Extension stops: 0°, 10°, 20°, 30° & 40°

Flexion stops: 45°, 60°, 75° & 90°

(brace is shipped with no flexion stops installed)

To change hinge stops, first remove the condyle pads held in place by hook/loop fasteners. Simply pull the pads off the hinge plates. To replace the pads be sure that they are positioned properly and firmly push them onto their hook/loop fasteners. Next, use a screwdriver to remove screws. Extension stops fit into the recessed slot on the face of the hinge plate, Flexion stops slide into the slot between the hinge plates from the back side. Remove existing stop, insert new stop, replace screw. Be sure that screws are seated and snug. Both medial (inside) hinge and lateral (outside) hinge stops must be identical in degrees.

### Care and cleaning

Straps and pads must be removed and cleaned regularly depending upon your activity level and perspiration. Hand wash straps and frame pads in cold water with mild detergent, rinse, AIR DRY. Occasionally apply a dry lubricant to hinge especially after use in water.

**B.** The frame pads are held in place by hook/loop fasteners. Carefully pull the pads off the brace frame. To replace the pads be sure that they are positioned properly and firmly push them onto their hook/loop fasteners.

The knee hinge condyle pads are held in place by hook/loop fasteners. Simply pull the pads off the hinge plates. To replace the pads be sure that they are positioned properly and firmly push them onto their hook/loop fasteners.

**C.** If you require a more snug fit of the brace at the knee, shims can be installed behind the condyle pads. Pull the condyle pads off the hinge plates, then place the shim with the loop side against the inner hinge plate, then replace the pad over the shim. Be sure the pads are positioned properly and firmly push them onto their hook/loop fasteners.

**D.** To remove the brace straps, loosen the hook/loop tabs at each strap end and feed the straps through the D rings on the medial (inside) side of the brace until the strap is free. Be sure to replace the straps according to their strap number locations, See illustration.

**E.** Strap 3, (the front strap) must be replaced as shown. Both ends of the strap will thread through the D rings from the inside and fold over onto the outside. The strap end with the rectangular piece of hook folds across the front side of the strap first. Then the strap end with the plastic tab closes over the front of the strap.

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