Latitude™ OA Knee Brace

INSTRUCTIONS FOR USE

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
©2023 Breg, Inc. All rights reserved.
Latitude is a trademark of Breg, Inc.

EC REP

E/U authorized representative
MDSS GmbH
Schiffgraben 41
D-30175 Hannover
Germany


Latitude OA Knee Brace

INSTRUCTIONS FOR USE



Table of Contents



Illustrations.....	3
<i>Languages:</i>	
English.....	7
Deutsch.....	8
Français.....	9
Español.....	10
Italiano.....	11
Norsk.....	12
Dansk.....	13
Svenska.....	14
Suomi.....	15
Nederlands.....	16
Português.....	17
日本語.....	18
简体中文.....	19
한국어.....	20
Polski.....	21



R = Right and L = Left
 R = Rechts und L = Links
 R = Droite et L = Gauche
 R = Derecha, L = Izquierda
 R = Destro e L = Sinistro
 R = Høyre og L = Venstre
 R = Højre og L = Venstre
 R = Höger och L = Vänster
 R = Oikea, L = Vasen
 R = Rechts, L = Links
 R = Direita e L = Esquerda
 R = 右, L = 左
 R = 右, L = 左
 R = 우측 및 L = 좌측
 R = strona prawa i L = strona lewa

Key
 Schlüssel
 Clé
 Llave
 Chiave
 Nøkkel
 Nøgle
 Nyckel
 Avain
 Sleutel
 Chave
 キー
 钥匙
 키
 Kluczyk

English

For single patient use only.

To be used under the supervision of a healthcare professional.

INITIAL FITTING PROCESS

1. Push the release button that says "PUSH" and slide to detach strap. Repeat for all four straps.
2. Sit on a stable surface with knee bent at approximately 45 degrees. Position brace on the lateral (outside) of the leg. Align the center of the hinge to slightly above the patella (knee cap).
3. Reattach the straps in sequence (1-4) by sliding clip into place.
4. For final brace adjustment, verify proper hinge alignment and adjust strap tension as needed.

MEDIAL / LATERAL OFFLOADING

5. **Medial offloading:** For medial compartment offloading, turn key counter-clockwise ("IN"). Drawing 5 shows a right leg with medial offloading.
6. **Lateral offloading:** For lateral offloading, turn key clockwise ("OUT"). Drawing 6 shows a right leg with lateral offloading.

BRACE COMFORT ASSESSMENT

Once the brace is applied, stand and walk a short distance to ensure comfort and proper correction. If there is any discomfort, adjust the straps as needed. If you experience pain in the knee joint or leg, turn the hinge adjustment back slightly. If more relief is required, turn the hinge adjustment no more than 1/4 of a turn, as directed by your healthcare professional.

⚠ WARNING: The brace is not designed to provide support for previous ligament injuries.

⚠ WARNING: If you experience increased pain, swelling, skin irritation, or any adverse reactions while using this product, immediately consult your medical professional.

⚠ WARNING: This device will not prevent or reduce all injuries. Proper rehabilitation and activity modification are also an essential part of a safe treatment program. Consult your licensed healthcare professional regarding safe and appropriate activity level while wearing this device.

CAUTION: Carefully read fitting instructions and warnings prior to use. To ensure proper performance of the brace, follow all instructions. Failure to properly position the brace and fasten each strap will compromise performance and comfort.

CAUTION: Use in high impact activities may cause premature wear of components.

CLEANING

Pads and straps may be hand washed in cold water and mild detergent. Rinse and lay flat to air dry. Do not machine wash or tumble dry. To clean the brace frame and hinge, wipe with a damp cloth. Avoid brace contact with salt water, sand, dirt and chlorinated water. Use a dry silicone spray if a hinge lubricant is needed.

If any serious incident has occurred in relation to the device, report to Breg and the competent authority of the Member State in which the user and/or patient is established.

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