ENGLISH

INTENDED USE/INDICATIONS:

The intended use of the cervical collars is to provide stabilization and motion restriction for the treatment of trauma-related injury, acute post-surgical recovery, and the management of chronic conditions of the cervical spine. Indicated for Post-surgical Stabilization, Decompression Procedures, Cervical Disc Syndrome, Fracture Management, Sprains / Strains, Cervical Stenosis, Cervicogenic Headaches, Trauma / Whiplash.

INTENDED USERS:

The intended user should be a licensed professional, the patient, or patient's caregiver. The user should be able to read, understand and be physically capable to perform all the directions, warnings and cautions provided in the information for use.

CONTRAINDICATIONS:

Not Applicable.

⚠️ WARNINGS AND PRECAUTIONS:

- Carefully read fitting instructions and warnings prior to use. To ensure proper performance of the brace, follow all instructions.
- If you experience, or have an increase in pain, swelling, skin irritation, or any adverse reactions while using this product, immediately consult your licensed healthcare professional.
- This device will not prevent or reduce all injuries. Proper rehabilitation and activity modification are also an essential part of a safe treatment program. Consult your licensed healthcare professional regarding safe and appropriate activity levels while wearing this device.
- Federal law restricts this device to sale by or on the order of a licensed healthcare practitioner.
- The frequency and duration of use shall be prescribed by a healthcare professional.
- For single patient use only.

INSTRUCTIONS FOR USE:

MEASUREMENT

1. Ensure the patient's head is in the desired position for accurate fitting.
2. Use a sizing guide to determine the correct collar size based on the trapezius muscle and chin alignment.

BACK PANEL POSITIONING

3. Bend the ends of the back panel to fit the neck shape.
4. Position the back panel under the neck with loop straps between the ears and the shoulders.

Note: Ensure the triangular cutout points down toward the patient's back.

FRONT COLLAR POSITIONING

5. Bend the ends of the front collar inward to conform to the neck.
6. Align the bottom of the collar with the patient's chest and attach straps to the sides.
7. Ensure equal tension on both straps.

TIP: Adjust side panels if chin piece presses on the throat.

USE AND CARE:

Hand wash in 30°C (86°F) water with mild soap, rinse thoroughly. Rinse and lay flat to air dry. Do not machine wash or tumble dry.

MATERIAL CONTENT:

Polypropylene, cotton, Hemp, High Density Polyethylene (HDPE), Nylon

WARRANTY:

Breg products come with a limited warranty against manufacturing and material defects. No other warranty, expressed or implied, is given and no affirmation of or by seller, by words or action, will constitute a warranty. Please see www.breg.com

DISPOSAL:

Dispose of product according to local regulations.

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