Breg.com/ascendpediatriccollar

SVENSKA

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ENGLISH

INTENDED USE/INDICATIONS:

The intended use of the Ascend pediatric collars is to provide stabilization and motion restriction for the treatment of trauma related injury, acute post-surgical recovery, and the management of chronic conditions of the cervical spine. Indicated for post-surgical stabilization, cervical disc syndrome, cervical stenosis, fracture management, sprains / strains, decompressive procedures, cervicogenic headaches, radiculopathy, trauma / whiplash.

INTENDED USERS:

The intended user should be a licensed professional, the patient, or patient's caregiver. The user should be able to read, understand and be physically capable to perform all the directions, warnings and cautions provided in the information for use.

CONTRAINDICATIONS:

Not Applicable.

⚠ WARNINGS AND PRECAUTIONS:

- Carefully read fitting instructions and warnings prior to use. To ensure proper performance of the collar, follow all instructions.
- If you experience, or have an increase in pain, swelling, skin irritation, or any adverse reactions while using this product, immediately consult your licensed health care professional.
- This device will not prevent or reduce all injuries. Proper rehabilitation and activity modification are also an essential part of a safe treatment program. Consult your licensed health care professional regarding safe and appropriate activity level while wearing this device.
- Federal law restricts this device to sale by or on the order of a licensed healthcare practitioner.
- Follow the licensing healthcare practitioner's protocol regarding method, frequency and duration of administration.

INSTRUCTIONS FOR USE:

BACK PANEL POSITIONING

- Bend the ends of the back panel inward to fit the neck shape.
- Position the patient's head in a neutral position and slide the back panel under the neck.
- Ensure loop straps are between the ears and shoulders, centering the back panel along the neck's mid line.

FRONT COLLAR POSITIONING

- Bend the ends of the front collar inward to conform to the neck.
- Place the collar against the patient's chest and align with the back panel straps.
- Attach straps securely, ensuring equal tension and proper alignment.
- TIP: Adjust side panels if chin piece presses on the throat.

USE AND CARE:

Hand wash in 30°C (86°F) water with mild soap, rinse thoroughly. Rinse and lay flat to air dry. Do not machine wash or tumble dry.

MATERIAL CONTENT:

Cotton, High Density Polyethylene (HDPE), Nylon

WARRANTY:

Breg products come with a limited warranty against manufacturing and material defects. No other warranty, expressed or implied, is given and no affirmation of or by seller, by words or action, will constitute a warranty. Please see www.breg.com

DISPOSAL:

Dispose of product according to local regulations.

NOTICE TO USERS:

If any serious incident has occurred in relation to the device, report to Breg and the competent authority of the Member State in which the user and/or patient is established.

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