ENGLISH

INTENDED USE/INDICATIONS:

The intended use of the Pinnacle Cervical Collar is to provide stabilization and motion restriction for the treatment of trauma related injury, acute post-surgical recovery, and the management of chronic conditions of the cervical spine. Indicated for post-surgical stabilization, cervical disc syndrome, cervical stenosis, fracture, managment, sprains / strains, decompressive procedures, cervicogenic headaches, radiculopathy, trauma / whiplash.

INTENDED USERS:

The user should be able to read, understand and be physically capable to perform all the directions, warnings and cautions provided in the information for use.

CONTRAINDICATIONS:

Not Applicable.

⚠ WARNINGS AND PRECAUTIONS:

- Carefully read fitting instructions and warnings prior to use. To ensure proper performance of the collar, follow all instructions.
- If you experience, or have an increase in pain, swelling, skin irritation, or any adverse reactions while using this product, immediately consult your licensed healthcare professional.
- This device will not prevent or reduce all injuries. Proper rehabilitation and activity modification are also an essential part of a safe treatment program. Consult your licensed healthcare professional regarding safe and appropriate activity level while wearing this device.
- This device is offered for sale by or on the order of a licensed healthcare professional. Follow application and care process described in instructions.
- For single patient use only.
- Federal law restricts this device to sale by or on the order of a licensed healthcare practitioner.
- Follow the licensing healthcare practitioner's protocol regarding method, frequency and duration of administration.

INSTRUCTIONS FOR USE:

BACK PANEL POSITIONING

1. Ensure patient's head is in a neutral position.
2. Slide back panel under patient's neck.
3. Ensure loop straps are positioned between patient's ears and shoulders. The back panel should be centered along the midline of the back of the neck.

FRONT COLLAR POSITIONING

4. Bend the back ends of the front collar inward to conform to the shape of the neck.
5. The front collar should start with the chin piece in its lowest position.
6. Place the bottom of front collar against the patient's chest and align sides with back panel loop straps.

CHIN ADJUSTMENT

7. While holding the front collar in position, pull out and rotate the knob clockwise to raise the chin piece until the patient's chin is supported and head is in a neutral position.
8. Once adjustments have been made, ensure the knob is secured into one of the stop positions.

STRAP ATTACHMENT

9. Attach straps to the hook on the sides of the front collar.
10. Ensure both straps have an equal amount of tension and proper alignment is maintained, resulting in a snug and supported fit.

NOTE: Ensure that chin piece is not pressing on patient's throat. If it is, adjust side panels of the front collar so they are off the shoulders.

USE AND CARE:

Hand wash in 30°C (86°F) water with mild soap, rinse thoroughly. Rinse and lay flat to air dry. Do not machine wash or tumble dry.

MATERIAL CONTENT:

Cotton, DMNT Resin, High Density Polyethylene (HDPE), Nylon

WARRANTY:

Breg products come with a limited warranty against manufacturing and material defects. No other warranty, expressed or implied, is given and no affirmation of or by seller, by words or action, will constitute a warranty. Please see www.breg.com

DISPOSAL:

Dispose of product according to local regulations.

NOTICE TO USERS: If any serious incident has occurred in relation to the device, report to Breg and the competent authority of the Member State in which the user and/or patient is established.

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