# EPIC™ LP LSO 637/650 LUMBAR ORTHOSIS

MODEL #SP30637-000

INCLUDES PADS

## INSTRUCTIONS FOR USE

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[www.breg.com](http://www.breg.com)



©2023 Breg, Inc. All rights reserved.  
Epic is a trademark of Breg, Inc.

AW-1.04450 Rev B 05/23

E/U authorized representative  
MDSS GmbH, Schiffgraben 41  
D-30175 Hannover, Germany

Size	1	2	3	4	EXT.
Waist Size	24" – 34" (61.0–86.4cm)	33" – 40" (83.8–101.6cm)	39" – 46" (99.1–116.8cm)	45" – 50" (114.3–127cm)	50" + (127cm+)

## ENGLISH

### MEASUREMENT

1. Measure patient's waist and refer to sizing chart (above) for corresponding size.

### ADJUSTING THE WAIST BELT SIZE

**Note:** The initial configuration is set for size 2.

2. To adjust the waist belt size, open the brace and lay it face down.
3. Flip the waist belt over to find the sizing numbers. Disengage the hook and loop and move one belt end to the desired sizing number.
4. Repeat on the other belt end to complete the sizing process.

### CLINICAL TIP

1. Flip waist belt over and apply to patient's waist.
2. Have patient hold the front of the waist belt to estimate appropriate length for sizing.
3. Fitter-adjust belt length to patient's size.
4. Flip waist belt back over and apply to patient.

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





## ENGLISH

### APPLYING THE BRACE

Before applying the brace, make sure the pull handles are positioned on the belt, close to the mesh cover on the back of the brace.

1. Align the center of the back pad to the patient's spine, and align the bottom edge of the back panel over lumbosacral joint. **Refer to picture 1 for proper orientation.**
2. Fasten the brace around the patient's waist by pulling the right belt over the left belt. Angle belts up or down to accommodate patient's anatomy.
3. Adjust the tension of the brace by pulling the pull handles away from the body.
4. Attach the pull handles to the belt.

### REMOVING THE BRACE

5. Detach the pull handles from the belt and position them close to the back of the brace.
6. Detach hook and loop at the front to remove the brace.

The brace can be trimmed, bent, molded, assembled, or otherwise customized to fit a specific patient. See product website for further instructions.

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 WARNINGS 

**WARNING:** Carefully read fitting instructions and warnings prior to use. To ensure proper performance of the brace, follow all instructions.

**WARNING:** If you experience, or have an increase in pain, swelling, skin irritation, or any adverse reactions while using this product, immediately consult your licensed health care professional.

**WARNING:** This device will not prevent or reduce all injuries. Proper rehabilitation and activity modification are also an essential part of a safe treatment program. Consult your licensed health care professional regarding safe and appropriate activity level while wearing this device.

**CAUTION:** This device is offered for sale by or on the order of a licensed health care professional.

**CAUTION:** Follow application and care process described in instructions.

**CAUTION:** For single patient use only.

**CARE & WASHING INSTRUCTIONS**

Hand wash in cold water with a mild detergent, rinse and lay flat to air dry.

**If any serious incident has occurred in relation to the device, report to Breg and the competent authority of the Member State in which the user and/or patient is established.**

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