
## POLAR CARE WAVE™

INSTRUCTIONS FOR USE  
INSTRUCCIONES DE USO  
MODE D'EMPLOI  
ISTRUZIONI PER L'USO  
GEBRAUCHSANWEISUNG  
GEBRUIKSAANWIJZING  
ANVÄNDARINSTRUKTIONER  
BRUKSANVISNING  
KÄYTTÖOHJEET  
BRUGSANVISNING  
INSTRUKCJA UŻYTKOWANIA  
ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ  
HASZNÁLATI UTASÍTÁS

**Breg, Inc.**  
2382 Faraday Avenue, Suite 300  
Carlsbad, CA 92008 USA  
Toll Free Tel: 1-800-321-0607  
Local: +1-760-795-5440  
Fax: +1-760-795-5295  
[www.breg.com](http://www.breg.com)

© 2024 Breg, Inc. All rights reserved.  
Polar Care Wave is a registered  
trademark of Breg, Inc.  
AW-1.04509 Rev C 05/24

E/U authorized  
representative  
MDSS GmbH  
Schiffgraben 41  
D-30175 Hannover  
Germany

UK Responsible Person  
Joint Operations Ltd  
11 Interface Business Park  
Royal Wootton Bassett  
Wiltshire  
England  
SN4 8SY

Scan QR Code for  
instructional video.

# POLAR CARE WAVE™

## INSTRUCTIONS FOR USE

**CAUTION:** United States Federal Law restricts this device to sale by or on the order of a medical professional.

The Polar Care Wave and cold compression pads are intended for single patient use only, dispensed with a prescription.

For questions or product inquiries, contact **Breg Customer Care** at **1-800-321-0607** or **+1-760-795-5440**.

### Table of Contents

Indications / Contraindications.....	2
Symbols Key.....	4
Polar Care Wave Description.....	5
Warnings.....	6
Cautions.....	9
Operating Instructions.....	10
Cleaning.....	15
Storage.....	16
Product Part Numbers / Warranty.....	17
Troubleshooting Guide.....	18
Product Specifications (Electrical, Environmental, Performance, and Physical).....	20
Electrical Safety.....	22
Discharge Reminders for Healthcare Professionals.....	23
Cold Compression Protocol.....	24

## INDICATIONS FOR USE

The Polar Care Wave System is intended to treat post-surgical and acute injuries to reduce edema, swelling, and pain where cold and compression are indicated. It is intended to be used by or on the order of licensed healthcare professionals in hospitals, outpatient clinics, athletic training settings, or home settings.

## CONTRAINDICATIONS

The patient should NOT use the Polar Care Wave System if the patient is suspected of or observed to have any of the following pre-existing conditions:

### Cold and Compression Therapy:

- |  |  |
|--|--|
| <ul style="list-style-type: none"> <li>• Significant vascular impairment in the affected region (e.g., from prior frostbite, arteriosclerosis, arterial insufficiency, diabetes, vascular dysregulation, or other vascular ischemic disease).</li> <li>• Current clinical signs of inflammatory phlebitis, venous ulcers, or cellulitis.</li> <li>• Any significant risk factors or current clinical signs of embolism (e.g., pulmonary embolus, pulmonary edema, cerebral infarction, atrial fibrillation, endocarditis, myocardial infarction, or atheromatous embolic plaque).</li> </ul> | <ul style="list-style-type: none"> <li>• A condition in which increased venous or lymphatic return is not desired in the affected extremity (e.g., lymphedema after breast cancer or other local carcinoma and/or carcinoma metastasis in the affected extremity).</li> <li>• Hypertension, cardiac failure, extreme low blood pressure, or decompensated cardiac insufficiency.</li> <li>• Localized unstable skin condition (e.g., dermatitis, vein ligation, gangrene, or recent skin graft) in the affected region.</li> <li>• Had recent toe surgery in the affected region.</li> </ul> |
|--|--|

### Compression Therapy:

- |  |  |
|--|--|
| <ul style="list-style-type: none"> <li>• Current clinical signs in the affected region of significant peripheral edema (e.g., deep vein thrombosis, chronic venous insufficiency, acute compartment syndrome, systemic venous hypertension, congestive heart failure, cirrhosis/liver failure, renal failure).</li> <li>• An acute, unstable (untreated) fracture in the affected region.</li> <li>• Any active local or systemic infection.</li> <li>• Obtunded or with diabetes mellitus, multiple sclerosis, poor circulation, spinal cord injuries, and rheumatoid arthritis.</li> </ul> | <ul style="list-style-type: none"> <li>• Areas of skin breakdown or damage (damaged or at-risk skin) producing uneven heat conduction across the skin (e.g., open wound, scar tissue, burn or skin graft). Any open wound must be dressed prior to use of the Polar Care Wave System.</li> <li>• Presumptive evidence of congestive heart failure</li> <li>• Pre-existing DVT condition</li> <li>• Deep acute venal thrombosis (Phlebothrombosis)</li> <li>• Inflammatory phlebitis process</li> </ul> |
|--|--|

Contraindications continued through page 4.


**CONTRAINDICATIONS****Compression Therapy (continued from page 2):**


- |  |  |
|--|--|
| <ul style="list-style-type: none"> <li>• Episodes of pulmonary embolism</li> <li>• Pulmonary edema</li> <li>• Acute inflammation of the veins (Thrombophlebitis)</li> <li>• Decompensated cardiac insufficiency</li> <li>• Arterial dysregulation</li> <li>• Erysipelas</li> <li>• Carcinoma and carcinoma metastasis in the affected extremity</li> </ul> | <ul style="list-style-type: none"> <li>• Decompensated hypertonia</li> <li>• Acute inflammatory skin diseases or infection</li> <li>• Venous or arterial occlusive disease</li> <li>• Medical situations where increased venous or lymphatic return is undesirable</li> <li>• Poor peripheral circulation</li> <li>• Severe arteriosclerosis, or active infection</li> </ul> |
|--|--|

**Cold Therapy:**

- |  |   |
|--|---|
| <ul style="list-style-type: none"> <li>• Known hematological dyscrasias that predispose to thrombosis (e.g., paroxysmal cold hemoglobinuria, cryoglobulinemia, sickle-cell disease, serum cold agglutinins).</li> <li>• Tissues inflamed as a result of recent injury or exacerbation of chronic inflammatory condition.</li> <li>• Compromised local circulation or neurologic impairment (including paralysis or localized compromise due to multiple surgical procedures or diabetes) in the affected region.</li> <li>• Cognition or communication impairments that prevent them from giving accurate and timely feedback.</li> <li>• Raynaud's disease or cold hypersensitivity (cold urticaria).</li> <li>• Decompensated hypertonia in the affected region.</li> <li>• Raynaud's phenomenon or other vasospastic conditions.</li> </ul> | <ul style="list-style-type: none"> <li>• Cold allergy</li> <li>• Cold agglutinin disorders like paroxysmal cold hemoglobinuria</li> <li>• Buerger's disease</li> <li>• Chilblains</li> <li>• Cryoglobulinemia</li> <li>• Sickle cell anemia</li> <li>• Diabetes</li> <li>• Hypersensitivity to cold</li> <li>• History of cold injury</li> <li>• Severe cardiovascular disease, anesthetic skin, hypercoagulation disorders, poor circulation, extremities sensitive to pain, extremely low blood pressure that are incapacitated, decreased skin sensitivity, vein ligation or recent skin grafts, or pheochromocytoma.</li> </ul> |
|--|---|

## SYMBOLS KEY

Label Symbols	Description
	This is the safety alert symbol. It is used to alert of potential personal injury hazards. Obey all safety messages that follow this symbol to avoid possible injury or death.
	Type BF, applied part
	Class II Equipment
	Refer to instructions for use
IP22	Ingress Protection Marking of IP22: Protection against solid foreign objects of 12.5mm and greater and against vertically falling water drops when tilted up to 15°.
	Disposal of product according to local regulations.
	Breg logo
	Lid Lock/Unlock
	Single Patient multiple use
	Lid to container alignment

Keypad Symbols	Description
	Colder setting
	Cold setting
	Regular compression setting
	Low compression setting
	On/off button


## POLAR CARE WAVE DESCRIPTION

The Polar Care Wave is a motorized cold therapy and intermittent compression device. Motorized cold therapy and compression have repeatedly been shown to decrease postoperative pain, swelling, inflammation, and narcotic use following a variety of surgical procedures, such as surgery to the shoulder, knee, and back.


The Polar Care Wave unit consists of a container, tubing, lid with pumps, power supply, and a cold compression pad. The cold compression pad is applied to the patient's injured area which delivers cold and intermittent compression therapy. The unit delivers cold therapy by circulating water from the container through the tubing into the cold compression pad. The compression therapy is delivered by the unit's air pump, which inflates and deflates air into the cold compression pad providing gentle squeezing of the pad around the affected area. The cold and compression therapies can be used in conjunction with each other or independently of each other, if preferred or prescribed by a licensed healthcare professional. Breg offers a family of cold compression pads to meet treatment needs (see page 12 for the full cold compression pad options).

The Polar Care Wave can be used and operated by a patient or caretaker in a home or clinical setting.


### These items are included:


### Accessories (sold separately):

Breg Ice Bags

Knee Y Adapter

Power Adapter

**⚠ WARNINGS****⚠ WARNING**

The Polar Care Wave can be cold enough to cause serious injury, including full skin necrosis. Follow and carefully read the instructions for use and the cold compression pad fitting instructions (printed on each cold compression pad) prior to use.

**General Warnings**

Cognitive risk factors	<p>Patients with any of the following cognitive risk factors should only use Polar Care Wave under direct supervision of a medical professional or caretaker. <b>If patient has any of the cognitive risk factors below, medical professional or caretaker should provide skin checks:</b></p> <ul style="list-style-type: none"> <li>• Children under 18 years old and elderly</li> <li>• Cognitive disabilities</li> <li>• Communication barriers</li> <li>• Use of medications that have a negative effect upon mental capacity</li> <li>• Visually impaired</li> </ul>
Discuss treatment with your licensed healthcare professional	<p>Ask your healthcare professional about potential adverse reactions and cold-induced injuries. Certain medical conditions make cold-induced injury more likely.</p>
Use only as prescribed	<p>Do not use this device if you did not receive or do not understand the instructions. Use only according to your healthcare professional's instructions regarding the frequency, duration, and settings of cold and compression application and length of breaks between uses, how and when to inspect the skin, and total length of treatment.</p> <p>United States Federal law restricts this device to sale by or on the order of a healthcare professional. The Polar Care Wave and cold compression pads are intended for single patient, multiple use only. Dispensed with a prescription. The prescription must include:</p> <ul style="list-style-type: none"> <li>• The frequency, duration, and settings of the cold and compression application</li> <li>• Length of breaks between uses</li> <li>• How and when to inspect the skin</li> <li>• Total length of treatment</li> </ul>

Warnings continued through page 9.

## **▲ WARNINGS**

<b>General Warnings</b>	
<b>Inspect skin regularly</b>	<p>Patients vary in sensitivity to cold, regularly inspect the skin under the cold compression pad (by lifting the edge) as prescribed, typically every 1 to 2 hours. Do not use the Polar Care Wave unit if skins checks are inhibited by a barrier.</p> <p>Stop using and consult your healthcare professional immediately if you experience any adverse reactions, such as: increased pain, burning, increased swelling, itching, blisters, increased redness, discoloration, welts, irritation and other changes in skin appearance, or any other reaction identified by your healthcare professional.</p>
<b>General safety</b>	<ul style="list-style-type: none"><li>• Keep tubing and cords away from the neck to avoid risk of strangulation.</li><li>• If a leak or a steady stream of water is flowing from any part of the unit, tubing or cold compression pad, discontinue use and contact Breg Customer Care. Excess moisture could result in a slip hazard or unwanted moisture at an injury site.</li><li>• Do not modify, service, or perform unspecified maintenance on this equipment or accessories. Only use parts and accessories designed for the Polar Care Wave.</li><li>• Do not attempt to sterilize Polar Care Wave or cold compression pads by any means.</li><li>• Set up and use this equipment in a location away from unsupervised children and pets.</li><li>• Do not operate the unit in a heavy dust/lint environment, excessive dust/lint may cause clogging. Do not leave the unit in direct sunlight, UV light may damage or discolor the unit.</li><li>• Device not intended for use with numbing agents.</li></ul>

### **NOTICE TO USERS**

If any serious incident has occurred in relation to the device, report to Breg and the competent authority of the Member State in which the user and/or patient is established.

**▲ WARNINGS****General Warnings**

Cold compression pad

- Regularly check the therapy area.
- Cold compression pads are non-sterile. Patients should take caution in applying cold compression pads over open wounds or breached skin. At a minimum, these areas should be cleaned and bandaged.
- Cold compression pad should be inspected for cleanliness for each treatment. If the pad is dirty, clean as indicated in the cleaning section on page 15.
- Check for moisture on the cold compression pad before placing on the skin. Remove any moisture before use.
- Ensure the therapy wrap is applied correctly before initiating any therapy. Do not wrap the cold compression pads as to restrict blood flow or fluid flow within the cold compression pad.
- Do not apply any casting or restrictive devices over the cold compression pad and tubing set as this may prevent proper operation and restrict skin checks.
- The Polar Care Wave and cold compression pads are intended for single patient, multiple use only. Re-use of single patient use pads may lead to risk of infection.

**Cold Therapy Warnings**

Cold risk

If the risk of cold-induced injury outweighs the benefits of cold therapy, do not prescribe Polar Care Wave.

Clinical risk factors

If the patient has any of the following clinical risk factors, use of Polar Care Wave may result in serious cold-induced injury, including full thickness skin necrosis:

- Pathologic sensitivity to cold
- Behaviors that negatively affect circulation, including poor nutritional status, smoking and tobacco use, excessive caffeine use, and excessive alcohol use
- Cold application area desensitization due to local anesthesia or regional nerve blocks
- Diabetes
- Hand/wrist or feet/ankle surgery
- Taking medications that have a negative effect on peripheral vascular circulation, including beta adrenergic blockers and local epinephrine use (such as in local anesthetics)

**▲ WARNINGS**

<b>Cold Therapy Warnings</b>	
<b>Cold risk mitigation</b>	<p>If you prescribe this product to patients with risk factors, consider taking special measures to control the risk, such as:</p> <ul style="list-style-type: none"> <li>• Recommending more frequent skin checks. If patient has cognitive risk factors, healthcare professional or caretaker should provide skin checks</li> <li>• Requiring more frequent follow-up examinations</li> <li>• Adding an insulation barrier between the cold compression pad and skin</li> <li>• Prescribing shorter durations of application, less frequent application, or eliminating nighttime application</li> <li>• Prescribe the warmer temperature setting (smaller snowflake/"cold" setting)</li> </ul>
<b>Application site</b>	Excessive moisture at the application site due to excessive bleeding, sweating, or condensation may increase the risk of serious cold-induced injury, including full thickness skin necrosis.




<b>Compression Therapy Warnings</b>	
<b>Compression risk</b>	If the risk of injury outweighs the benefits of compression therapy, do not prescribe Polar Care Wave.
<b>Clinical risk factors</b>	<p>Patients with any of the following compression risk factors should consult with a healthcare professional prior to use:</p> <ul style="list-style-type: none"> <li>• Neuropathies or tissue viability problems (i.e. diabetes, arterial or venous insufficiencies.)</li> <li>• Patients positioned in the supine lithotomy position for extended lengths</li> <li>• Poor circulation</li> </ul>
<b>Compression risk mitigation</b>	Immediately stop compression therapy if you experience any sense of discomfort, numbness or tingling of the limb and consult your healthcare professional.

**CAUTIONS**

- Do not walk with tubing connected to the cold compression pads. It is a trip hazard.
- If the unit is damaged or cracked, discontinue use and contact Breg Customer Care.
- Do not use hot water in this unit. It may be unsafe or damage the unit since it has not been designed or tested with use of hot water.
- Do not use any other fluids besides water or recommended cleaning agents on page 15 with the Polar Care Wave. Doing so will damage the unit.
- Do not run cold therapy without water. Using cold therapy without water can cause permanent damage to the water pump.



## OPERATING INSTRUCTIONS


## Setup with Water and Ice:

1	Set up and operate unit on a flat surface so the unit does not tip over.	
2	Unlock handle by rotating handle to the front of the unit.	
3	<p>Remove lid and fill container with water to indicated water fill line.</p> <p>Fill to lines on label inside container. Use cubed ice for optimal performance.</p> <p>Only use ice and water or four (4) Breg Ice Bags and water with this unit (see page 11 for Setup with Breg Ice Bag Accessories).</p> <p><b>CAUTION:</b> Do not overfill the unit with ice and water. Appropriate levels are indicated by the ice/water label in the container. Overfilling the unit can cause water to spill out or prevent the lid from closing properly.</p>	
4	Place lid on container using align stickers.	
5	Lock by rotating the handle to the upright or back position according to image on handle.	
6	Every 6-8 hours or as needed, pour out enough water to return water level to original water fill line and re-fill ice to ice line.	

Continue to Power Supply instructions on page 11.






## OPERATING INSTRUCTIONS

Setup with Breg Ice Bag Accessories and Water:	
<p><b>1</b> Remove eight (8) empty Breg Ice Bags from the packaging. Only use Breg Ice Bags that come with the Polar Care Wave.</p>	
<p><b>2</b> Fill the four eight (8) Ice Bags with water to fill line. Do not over fill.</p>	
<p><b>3</b> Freeze Breg Ice Bags in freezer. Dry Breg Ice Bag exterior before placing in freezer.</p>	
<p><b>4</b> Set up and operate unit on a flat surface so the unit does not tip over.</p>	
<p><b>5</b> Unlock handle by rotating handle to the front of the unit. Remove lid and add water to fill line in container.</p>	<div style="border: 1px solid black; padding: 5px; text-align: center;"> <p><b>FIRST ADD WATER SECOND ADD ICE BAGS</b></p> </div> 
<p><b>6</b> Place four (4) Breg Ice Bags into the corners of the container as shown. Only use ice and water or four (4) Breg Ice Bags and water with this unit.</p>	
<p><b>7</b> Replace lid and lock by rotating the handle to the upright or back position.</p>	
<p><b>8</b> Replace four (4) melted Breg Ice Bags with four (4) frozen Breg Ice Bags every 4-5 hours or as needed. While using the Breg Ice Bags, cubed ice may also be added to the container up to the ice fill line to optimize the colder setting.</p>	

Power Supply:	
<p><b>1</b> Plug power supply into the wall and other end of cord into lid.</p>	


## OPERATING INSTRUCTIONS

**Cold compression Pad Application:****(Make sure Breg logo and pad instructions are facing out, away from the injury site)**


Standard / Large Knee Pad	
Hip Pad	
Shoulder Pad	
Back Pad	
Foot / Ankle Pad	
Universal Pad (use for knee or shoulder)	

## OPERATING INSTRUCTIONS

## Connecting and Disconnecting Cold compression Pad to Unit:


<p><b>1</b> Line up the Breg logos on the tubing and pad connectors.</p>	
<p><b>2</b> Squeeze down gently on the blue clips.</p>	
<p><b>3</b> Push connectors together until they make a "click" sound.</p>	
<p>Disconnect the cold compression pad from the unit by pressing the blue clips and gently pulling the connectors apart.</p> <p><b>4</b> The Breg Polar Care cold compression pad may be disconnected from the unit while it's running and connectors will seal, preventing water from flowing out of tubing. <b>Note:</b> Some dripping during release is normal.</p>	

## Power Adapter:

	<p>First, take the power supply and appropriate plug adapter from the packaging.</p> <p>Then with the power supply plugs folded in vertically (shown on left), align the power supply and plug adapter (shown on right). Make sure locking notch is on top, otherwise invert the plug adapter.</p> <p>Finally, slide the adapter down gently until it is fully recessed down to the power supply.</p>
<div style="display: flex; justify-content: space-around; align-items: center;"> <div style="text-align: center;">  <p><b>TYPE C</b></p> </div> <div style="text-align: center;">  <p><b>TYPE I</b></p> </div> <div style="text-align: center;">  <p><b>TYPE G</b></p> </div> </div>	

## OPERATING INSTRUCTIONS

### Keypad Key:

### Operating the Unit:

Press the on/off button to turn the unit on. It will illuminate.

**1**

Upon initial start-up, the unit will default to colder and regular settings for cold and compression therapy, respectively.

If the unit was previously in use, it will default to the last settings selected when the unit was turned off.

**2**

Press the cold button to cycle between colder, cold, and off.

**3**

Press the compression button to cycle between regular, low, and off.

**4**

Press the on/off button to turn the unit off. All lights will turn off.

- Cold and compression therapies may be used at the same time or independently.
- If the on/off button continuously flashes, refer to troubleshooting guide on page 18.
- If the on/off button is on and both therapies are off for more than 15 minutes, the unit will automatically shut off.
- Some condensation on the tube and cold compression pad is possible, especially in warmer climates.

## CLEANING

**▲ WARNING:** Do not apply cold compression pad to patient during the cleaning procedure. Exposure to cleaning agent may cause injury.

**CAUTION:** Do not apply a direct stream of any liquid onto the unit's lid, submerge the unit, or allow any liquid to pool on the top of the lid. The unit's lid contains all of the unit's electronics.

Clean as necessary or if soiled or debris is present in the unit. If using well water, more frequent cleanings may be required.

### Polar Care Wave – Exterior Surfaces (container, lid, tubing, cold compression pad, and y-adapter)



Cleaning agents for home setting:	<ul style="list-style-type: none"> <li>Mild soap</li> <li>Water</li> </ul>
Cleaning agents for clinical setting:	<ul style="list-style-type: none"> <li>Use only mild soap and water on cold compression pad</li> <li>Disposable disinfection wipes</li> <li>Quaternary ammonium</li> </ul>
Procedure:	<ul style="list-style-type: none"> <li>Follow the manufacturer's instructions and precautions for the cleaning agent you select.</li> <li>Use selected cleaner and a soft cloth to wipe surfaces of container, lid, tubing, connectors, and cold compression pad.</li> <li>Wipe dry.</li> </ul>

### Polar Care Wave – Internal Components (pump system, internal tubing, inside of cold compression pad and y-adapter)

Cleaning agents:	<ul style="list-style-type: none"> <li>Bleach</li> <li>Water</li> </ul>
Procedure:	<ul style="list-style-type: none"> <li>Connect cold compression pad to unit.</li> <li>Fill the unit with cold tap water to the ice line and add 1/3 cup of bleach (5 tablespoons or 79 milliliters).</li> <li>Do not add ice.</li> <li>Place the lid assembly on the cooler and run the cold therapy for at least 10 min.</li> <li>Stop the unit and turn off, disconnect the power supply and place the lid assembly in the sink or drain.</li> <li>Raise the pad above the unit and roll the cold compression pad while still connected and squeezing the water out of the pad, into the drain or sink.</li> <li>Disconnect the cold compression pad, pour the water out of the cooler, rinse and wipe dry.</li> </ul>

## STORAGE

**Storing the Unit:**

<b>1</b>	Ensure the unit is off and the power supply is disconnected.	
<b>2</b>	Remove the lid assembly from the unit and place it in the sink or drain.	
<b>3</b>	Remove water from the cold compression pad by rolling the pad up while still connected and squeezing the water out of the pad, into the drain or sink.	
<b>4</b>	Disconnect the cold compression pad.	
<b>5</b>	Drain the cooler and wipe dry the unit.	
<b>6</b>	Coil tubing as shown and secure using straps.	
<b>7</b>	Do not store with water. Refer to page 20 for environmental storage requirements.	

[breg.com/wave](http://breg.com/wave)

Scan QR Code for instructional video.

---

**PRODUCT PART NUMBERS**

---

Product #	Description
100716-000	PC Wave INTL
C00003	Breg Cold Compression Pad, Knee L
C00004	Breg Cold Compression Pad, Shoulder
C00005	Breg Cold Compression Pad, Foot/Ankle
C00013	Breg Cold Compression Pad, Hip
C00016	Breg Cold Compression Pad, Universal
C00017	Breg Cold Compression Pad, Knee
C00020	Breg Cold Compression Pad, Back
C00015	Breg Carrying Bag
100576-000	PC Wave Knee Y Adapter
100582-000	Breg Ice Bags
100717-000	International PC Wave Replacement Power Supply
100718-000	Power Supply Adapter Kit, INTL

---

**WARRANTY**

---

Breg, Inc. warrants that this product is free from defects in workmanship and materials, and fit for use for 180 days from initial purchase under normal use for which it was intended.

Breg, Inc.'s obligation under this warranty is limited to the replacement or repair of any defective part or parts of this product. All expressed or implied warranties, including the warranty of merchantability and fitness for a particular purpose, are limited to the actual warranty period set forth above. No other warranty, expressed or implied, is given and no affirmation of or by seller, by words or action, will constitute a warranty.

## TROUBLESHOOTING GUIDE

Problem	Possible Solutions
On/off button flashing notification	<ul style="list-style-type: none"> <li>• Restart the Unit               <ul style="list-style-type: none"> <li>○ Turn unit off by either pressing the on/off button or by disconnecting and reconnecting power supply from the unit or from the wall.</li> <li>○ Press on/off button to turn unit back on.</li> <li>○ Press cold or compression therapy buttons (if required) to start desired therapies.</li> </ul> </li> <li>• If compression therapy does not function after pressing compression button (cold compression pad doesn't appear to be inflating or compression therapy setting light is not on), see Compression Therapy troubleshooting sections on page 19.</li> <li>• If cold therapy does not function after pressing cold button, see Cold Therapy troubleshooting section on page 19.</li> </ul>
On/off button isn't responding	<ul style="list-style-type: none"> <li>• Confirm power outlet is working and power supply is fully inserted in the outlet.</li> <li>• Confirm that power supply plug is fully inserted into the unit lid.</li> <li>• Confirm the correct Polar Care Wave power supply is connected (see page 20).</li> </ul>
Condensation	<ul style="list-style-type: none"> <li>• Wrap material loosely over cold compression pad and tubing to minimize air exposure.</li> <li>• Protect the wound site by using a sterile dressing with waterproof barrier. Acceptable barriers include an undercast bandage similar to Webril™ (ARTG 371503), a gauze bandage similar to Kerlix™ (ARTG 186369) or a tubular bandage.</li> </ul> <p><b>Note:</b> Some condensation on the tube and cold compression pad is possible, especially in warmer climates. Place a cloth layer under the container to protect items under the unit.</p>
Unit is leaking from blue connectors	<ul style="list-style-type: none"> <li>• Disconnect blue tabs on the cold compression pad connector.</li> <li>• Check that the Breg logos on the blue connectors are lined up correctly and reconnect.</li> </ul>

**If the problem continues, please contact Breg Customer Care at 1-800-321-0607 or +1-760-795-5440 for assistance.**

## TROUBLESHOOTING GUIDE

Problem	Possible Solutions
Cold therapy: water not flowing to cold compression pad or pad not cold	<ul style="list-style-type: none"> <li>• Confirm on/off button is illuminated and a cold therapy "colder" or "cold" setting light is on to indicate that cold therapy is functioning.</li> <li>• Check that blue connectors between cold compression pad and tubing are fully engaged with Breg logos lined up. Disconnect blue connectors and reconnect if needed.</li> <li>• Allow 10 minutes for cold compression pad temperature to stabilize.</li> <li>• Confirm unit has both ice and water to appropriate levels.</li> <li>• While using the Breg Ice Bags, cubed ice may also be added to the container up to the ice fill line to optimize the colder setting.</li> <li>• Elevating the unit to the same or higher elevation than the cold compression pad may increase water flow to the cold compression pad.</li> <li>• Check that tubing and cold compression pad are not kinked.</li> <li>• Remove cold compression pad and lay it flat. Allow cold compression pad to fill with water and then reapply.</li> </ul>
Compression therapy: cold compression pad not compressing or filling with air	<p><b>Note:</b> During compression therapy operation, it is normal for the compression to inflate and deflate continuously. The duration of each inflation-deflation cycle is approximately 1 to 2 minutes.</p> <ul style="list-style-type: none"> <li>• Confirm on/off button is illuminated and a compression therapy "regular" or "low" setting light is on to indicate that compression therapy is functioning.</li> <li>• Check that blue connectors between cold compression pad and tubing are fully engaged with Breg logos lined up. Disconnect blue connectors and reconnect if needed.</li> <li>• Tighten cold compression pad.</li> <li>• Check that tubing and cold compression pad are not kinked.</li> <li>• Remove cold compression pad and lay it flat. Check if the cold compression pad is inflating with air and then reapply.</li> </ul>
Compression therapy: cold compression pad not deflating properly	<ul style="list-style-type: none"> <li>• Check that blue connectors between cold compression pad and tubing are fully engaged with Breg logos lined up.</li> <li>• Disconnect blue connectors, deflate pad by rolling to remove air, and reconnecting blue connectors if needed.</li> </ul>

**If the problem continues, please contact Breg Customer Care at 1-800-321-0607 or +1-760-795-5440 for assistance.**

## PRODUCT SPECIFICATIONS

### Electrical Specifications (Also see Electromagnetic Immunity Test Table on page 21)

Unit Rated: 12VDC, 1A

Power supply is specified as part of ME Equipment

Manufacturer: HDP

Model: HDP12-MD12010C

Output: 12VDC, max. 1.0A

Back up supply: FRIWO

Model: FW8000M/12

Output: 12VDC, max. 1.0A

Rated Input: 100-240VAC, 50-60Hz

Standards Compliance:

Designed to conform to applicable requirements of: EN 60601-1, Edition 3.1 and 60601-1-2, 4th Edition

Electrical Classifications:

Class II, Continuous Operation. Not suitable with flammable anesthetics. Group 1, Class B.

### Environmental Requirements

Transport & Storage:

-25°C to 5°C, and

5°C to 35°C at a relative humidity of 10% to 90%, non-condensing;

>35°C to 70°C at a water vapor pressure up to 50 hPa

Operating:

5°C to 40°C (41°F to 104°F) 15-90% RH non-condensing,  
700 hPa to 1060 hPa

Indoor use only up to 2000m altitude

There is no warm up time required from when the Wave is stored at the minimum storage temperature and when the product is used.

### Performance Specifications

Expected Service Life:

1 year from date of first use (unit and accessories).

Average operating temperature of a Polar  
Care Wave cold compression pad:

≥45°F at the colder setting

≥50°F at the cold setting

**Note:** The Foot/Ankle pad at the cold setting may operate at temperatures above 60°F which only provides a cooling sensation.

Average peak compression pressure of the  
cold compression pad:

50 mmHg at the regular setting

25 mmHg at the low setting

**PRODUCT SPECIFICATIONS**
**Physical Specifications**

Weight: 5 lbs (2.3 kg) Empty, 16 lbs (7.3 kg) Filled

**Electromagnetic Immunity Test Table**

Immunity Test	Basic EMC Standard	IEC 60601 Test Level (Home Healthcare)
Electrostatic discharge (ESD)	IEC 61000-4-2	+/- 8 kV contact +/- 2kV, +/-4kV, +/-8 kV, +/- 15kV air
Radiated RF EM fields	IEC 61000-4-3	10V/m 80 MHz - 2.7 GHz 80% AM at 1 kHz
Proximity fields from RF wireless communications equipment	IEC 61000-4-3	Reference IEC 60601-1-2 table 9
Rated power frequency magnetic fields	IEC 61000-4-8	30 A/m 50Hz or 60Hz
Electrical fast transients/bursts	IEC 61000-4-4	+/- 2 kV for power supply lines 100 kHz repetition frequency
Surge	IEC 61000-4-5	+/- 0.5 kV, +/- 1 kV line-to-line +/- 0.5 kV, +/- 1 kV, +/- 2 kV line-to-ground
Conducted disturbances induced by RF fields	IEC 61000-4-6	3V 0.15MHz - 80 MHz 6V in ISM and amateur radio bands between 0.15 MHz and 80 MHz 80% AM at 1 kHz
Voltage dips	IEC 61000-4-11	0% UT 0.5 cycle At 0, 45, 90, 135, 180, 225, 270, and 315 degrees  0% UT: 1 cycle 70% UT: 25/30 cycles Single phase: at 0 degrees
Voltage interruptions	IEC 61000-4-11	0% UT: 250/300 cycles

**Note:** UT is the mains voltage prior to application of the test level

**⚠ ELECTRICAL SAFETY**

Flammables	<ul style="list-style-type: none"> <li>When using or cleaning the Polar Care Wave, avoid flammables and oxidizers. Do not use in places with flammable vapors or gases (e.g., flammable anesthetics), high oxygen concentrations, or other oxidizers (e.g., nitrous oxide).</li> </ul>
Keep electrical components dry and free of debris	<ul style="list-style-type: none"> <li>Even though the Polar Care unit is designed to hold and pump water, the electrical components of the unit and power supply must be kept dry.</li> <li>Do not handle the power supply or electrical cord with wet hands.</li> <li>Always keep the unit in a place where the connections will not fall into water (e.g., a tub, sink, etc.).</li> <li>If the electrical connections fall into water, do not touch any part of the product.</li> <li>Unplug only at dry electrical connections.</li> </ul>
Power supply	<ul style="list-style-type: none"> <li>Use only the power supply provided with the Breg Polar Care Wave unit. Model number is listed in Electrical Specification section on page 20.</li> <li>Keep the cord and power supply away from heated surfaces.</li> <li>Never operate this product if it has a damaged cord, power supply, or plug, which could result in electric shock.</li> <li>If the cord, power supply, or plug is damaged, unplug and contact Breg Customer Care.</li> </ul>
Electromagnetic interference and emissions	<ul style="list-style-type: none"> <li>This Polar Care Wave unit may cause electromagnetic or other interference with other electrical devices. Avoid using the Polar Care Wave unit adjacent to or stacked on other equipment because it could result in improper operation of the unit or the other equipment.</li> <li>Portable RF communications equipment should be used no closer than 12 inches (30 cm) to any part of the Polar Care Wave unit, doing so may result in product performance degradation.</li> <li>Using accessories or a power supply other than those provided could result in increased electromagnetic emissions or decreased electromagnetic immunity of the unit and can result in unit malfunction.</li> </ul>

**Disposal**

This product contains electrical components and can not be discarded in standard waste. Dispose of product according to local regulations.

## DISCHARGE REMINDERS FOR HEALTHCARE PROFESSIONALS

Follow these reminders prior to discharging the patient from facility care to home use.

- 1. Patient Screen.** Screen the patient for any contraindications and/or applicable warnings. If the patient has any contraindications (see Contraindications), do not dispense Breg's Polar Care Wave to the patient. If any of the warnings apply to the patient (see Warnings), determine the appropriateness of application of Breg's Polar Care to that patient.
- 2. Instructions For Use.** Instruct the patient on how to properly use Breg's Polar Care Wave. Review the Operating Instructions in this document and affixed to the unit with each patient.
- 3. Prescription.** Instruct the patient regarding the healthcare professional's prescribed protocol: frequency, compression level, temperature level and duration of use, length of breaks between uses, how and when to inspect the skin, and total length of treatment. A protocol template is included in this document (see page 24) to fill in and provide to the patient. The duration of application, temperature and compression level may vary depending upon the patient. If the patient does not experience pain relief, the physician may increase the duration of application and/or adjust temperature and compression levels. As the application duration is increased, the frequency of the skin inspections should increase (see Warnings). Instruct the patient to take breaks by turning off the unit for a specified period of time.
- 4. Potential For Injury.** Inform the patient that improper use can result in serious skin injury, including full thickness skin necrosis. Emphasize the importance of following the prescribed protocol, proper cold compression pad application, and skin inspection.
- 5. Proper Cold compression Pad Application.** Instruct the patient to use only the Breg Polar Care Wave cold compression pad designed for the body part being treated; other pads may not connect or may be colder, increasing the risk of serious cold-induced injury, including full thickness skin necrosis. Do not use the Polar Care Wave unit if skin checks are inhibited by a barrier.
- 6. Skin Inspection.** Instruct the patient to inspect the skin receiving treatment per the practitioner's instructions, typically every 1 to 2 hours. Do not use the Polar Care Wave unit if skin checks are inhibited by a barrier.
- 7. Discontinue.** Instruct the patient to stop using Breg Polar Care Wave and consult his/her healthcare professional immediately if he/she experiences any adverse reactions, such as: increased pain, burning, increased swelling, itching, blisters, increased redness, discoloration, welts, or other changes in skin appearance.
- 8. Documentation.** Provide the patient a prescription. The prescription must include:
  - The frequency, duration, and settings of the cold and compression application
  - Length of breaks between uses
  - How and when to inspect the skin
  - Total length of treatment

## Cold Compression Protocol\*

Treatment Period	Cold Frequency/Duration/Setting(s)	Compression Frequency/Duration/Setting(s)	Skin Inspection Frequency/Breaks Between Use
Week 1			
Weeks 2-4			
Week 4+			

\*To be completed by a licensed healthcare professional.

### Discuss treatment with your licensed healthcare professional

Ask your healthcare professional about potential adverse reactions and cold-induced injuries. Certain medical conditions make cold-induced injury more likely.

### Use only as prescribed

Do not use this device if you did not receive or do not understand the instructions. Use only according to your healthcare professional's instructions regarding the frequency, duration, and settings of cold and compression application and length of breaks between uses, how and when to inspect the skin, and total length of treatment.

United States Federal law restricts this device to sale by or on the order of a healthcare professional. The Polar Care Wave and cold compression pads are intended for single patient, multiple use only. Dispensed with a prescription. The prescription must include:

- The frequency, duration, and settings of the cold and compression application
- Length of breaks between uses
- How and when to inspect the skin
- Total length of treatment

### Inspect skin regularly

Patients vary in sensitivity to cold, regularly inspect the skin under the cold compression pad (by lifting the edge) as prescribed, typically every 1 to 2 hours. Do not use the Polar Care Wave unit if skins checks are inhibited by a barrier.

Stop using and consult your healthcare professional immediately if you experience any adverse reactions, such as: increased pain, burning, increased swelling, itching, blisters, increased redness, discoloration, welts, irritation and other changes in skin appearance, or any other reaction identified by your healthcare professional.

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