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

ENGLISH
INTENDED USE/INDICATIONS: Soft knee bracing products are designed to provide control and support for ligament injuries and instabilities. Indicated to support the knee and provide protected range of motion following mild to moderate MCL and/or LCL instabilities, mild to moderate ACL and/or PCL instabilities and meniscal repairs.
INTENDED USERS: The intended user should be a licensed professional, the patient or patient's caregiver. The user should be able to read, understand and be physically capable to perform all the directions, warnings and cautions provided in the information for use.
CONTRAINDICATIONS: Not Applicable.
WARNINGS AND PRECAUTIONS:
• If you experience increased pain, swelling, skin irritation, or any adverse reactions while using this product, immediately consult your medical professional.
• This device will not prevent or reduce all injuries. Proper rehabilitation and activity modification are also an essential part of a safe treatment program. Consult your licensed health care professional regarding safe and appropriate activity level while wearing this device.
• The United States federal law restricts this device to sale by or on the order of a licensed health care practitioner.
• For single patient use only.
• Carefully read fitting instructions and warnings prior to use. To ensure proper performance of the brace, follow all instructions. Failure to properly position the brace and fasten each strap will compromise performance and comfort.
• Use in high-impact activities may cause premature wear of components.
• Care and cleaning of this product is essential to its continuing strength and performance. Follow removal and cleaning process described in instructions.

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RANGE OF MOTION (ROM) STOP
CrossRunner Adjustment Instructions – Polycentric hinge with Stops:
Dual function extension and flexion ROM Stops
• Every Stop may be used to control Flexion or Extension
• 7 Unique Stops provide 14 angle set positions:
o 80 Flexion / 10 Extension (Pre-Installed)* o 70 Flexion / 20 Extension
o 45 Flexion / 45 Extension o 75 Flexion / 15 Extension
o 50 Flexion / 40 Extension o 90 Flexion / 0 Extension (Quantity 2)
o 60 Flexion / 30 Extension (Quantity 2)

ROM STOP REMOVAL

INSERTION OF THE ROM STOPS
• Select the ROM Stop
o Extension – Anterior side
o Flexion – Posterior side
• Press firmly to snap in ROM Stop

USE AND CARE:

• Attach Strap ends to strap Prior to cleaning.
• After exposure to salt water or dirt, thoroughly rinse the brace with fresh water and allow to air dry.
• Hand wash in 30°C (86°F) water with mild soap, rinse thoroughly. Rinse and lay flat to air dry. Do not machine wash or tumble dry.
• A dry lubricant, such as Teflon spray, may be used on the hinge occasionally. Do not use wet lubricants.
• The hinge assemblies may be removed to allow easier washing of the sleeve and straps.

MATERIAL CONTENT:

Acetal, Aluminum, BreathFi™, Nylon, Stainless Steel, Polyurethane
Breg products come with a limited warranty against manufacturing and material defects. No other warranty, expressed or implied, is given and no affirmation of or by seller, by words or action, will constitute a warranty. Please see www.breg.com

DISPOSAL:

Dispose of product according to local regulations.

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NOTICE TO USERS: If any serious incident has occurred in relation to the device, report to Breg and the competent authority of the Member State in which the user and/or patient is established.

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