
DVT Guardian™

Instructions For Use Instrucciones de uso Mode d'emploi

The product is composed of the control unit and the leg wrap cuff. The control unit can be charged by the power charger supplied. The control unit and the leg wrap cuff cannot be separated.

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www.breg.com/products/dvt-prophylaxis

2382 Faraday Avenue, Suite 300
Carlsbad, CA 92008 USA
Toll Free Tel: 1-800-321-0607
Local: +1-760-795-5440
Fax: +1-760-795-5295
www.breg.com

©2024 Breg, Inc. All rights reserved.
DVT Guardian is a trademark of Breg, Inc

11.000024 Rev B 1124

Rx Only

Table of Contents

Symbol Key.....	2
Introduction.....	3
Indications for Use.....	3
Contraindications.....	3
Warnings and Precautions.....	3
Technical Information.....	4
Condition of Working, Storage and Transportation.....	4
Operating Instructions.....	5
Cleaning and Disinfecting.....	7
Troubleshooting.....	7
Maintenance and Disposal.....	7
Warranty.....	7
Electromagnetic Compatibility.....	7

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Label Symbols	Description
	This is the safety alert symbol. It is used to alert of potential personal injury hazards. Obey all safety messages that follow this symbol to avoid possible injury or death.
	Type BF, applied part
	Class II Equipment
	Refer to instructions for use
IP22	Ingress Protection Marking of IP22: Protection against solid foreign objects of 12.5mm and greater and against vertically falling water drops when tilted up to 15°.
	Disposal of product according to local regulations.
	Temperature limit
	Single Patient multiple use
	Medical device symbol
	No presence of natural rubber latex
	ON/OFF button
Rx Only	This device is restricted for sale on or by the order of a physician

Introduction

The DVT Guardian is intended to be an easy-to-use portable system, prescribed by a physician, for use in the home by patients or clinical setting by clinicians to help prevent the onset of DVT in patients by stimulating blood flow in the extremities (stimulating muscle contractions).

The purpose of the device is to aid in the prevention of Deep Vein Thrombosis (DVT) by stimulating blood flow in the legs. This is accomplished by an electronically controlled pump delivering a preset amount of air to the leg wrap cuffs that, in turn, compress the calf to aid blood flow out of the lower extremities.

Indications for Use

The device is intended to be an easy to use portable system, prescribed by a physician, for use in the home or clinical setting to help prevent the onset of DVT in patients by stimulating blood flow in the extremities (stimulating muscle contractions). This device can be used to:

- Aid in the prevention of DVT;
- Enhance blood circulation;
- Diminish post-operative pain and swelling;
- Reduce wound healing time;
- Aid in the treatment and healing of stasis dermatitis, venous stasis ulcers, arterial and diabetic leg ulcers, chronic venous insufficiency and reduction of edema in the lower limbs.

The unit can also be used as an aid in the prophylaxis for DVT by persons expecting to be stationary for long periods of time.

⚠ Contraindications

The device **MUST NOT** be used to treat the following conditions:

- Persons with suspected, active or untreated: deep vein thrombosis, ischemic vascular disease, severe arteriosclerosis, pulmonary edema, severe congestive heart failure, thrombophlebitis, or an active infection.
- On the legs where cuffs would interfere with the following conditions: vein ligation, gangrene, dermatitis, open wounds, a recent skin graft, massive edema or extreme deformity of the leg.
- On any neuropathy.
- On extremities that are insensitive to pain.
- Where increased venous or lymphatic return is undesirable.
- Acute inflammatory skin diseases.
- Arrhythmia
- Erysipelas
- Unstable hypertension

⚠ Warnings and Precautions

- The device wraps are designed for single patient use only.
- Operation of the device can be done by the patient.
- To prevent extremity compartment syndrome, special attention should be given to patients who are positioned in the supine lithotomy position for extended lengths of time. This includes patients with or without wraps.
- Patients with diabetes or vascular disease require frequent skin assessment. Consult a Physician.
- The device used alone or in combination with warming device may cause skin irritation. Regularly check for discomfort, compliance, and skin irritation.
- Stop using the device if swelling occurs; consult a Physician.

- Only the original accessories can be used, especially the power charger. Use of other accessories may cause damage to the device and/or patient.
- Do not open or remove covers. No user serviceable parts inside. Direct all unit issues to the manufacturer or seller.
- Do not attempt to service the device while device is in use.
- If pulsations or throbbing occur, the wrap cuff may be wrapped too tightly. Loosen immediately.
- Device is to be used only by the patient prescribed, and only for its intended use.
- Ensure the device is turned off and unplugged from the wall outlet prior to and while cleaning or disinfecting.
- When charging, place the device in a position where it is easy to disconnect the power charger from the power network.
- Equipment should not be used in the presence of any flammable anesthetic mixture with air, oxygen, or nitrous oxide.
- Do not operate device in a wet environment.
- Do not immerse in any liquid for any reason.
- Do not place the device in an autoclave.
- Do not expose the device to the following conditions which may damage the device: direct sunlight, excessive lint or dust.
- Do not use or store the device outside the temperature range and humidity specified by the manufacturer, which may damage the product.
- Allow wraps to warm to room temperature if exposed to temperatures below 5 °C (41°F).
- Do not subject the device to extreme shocks, such as dropping the device.
- Keep the device away from infants and children. Pets and pests may damage the device. Check the device each time before use.
- Medical Electrical Equipment needs special precautions regarding Electromagnetic Compatibility (EMC). Portable and mobile Radio frequency (RF) Communications Equipment can be affected by other Medical Electrical Devices.
- Do not run the product under the environment of strong electromagnetic interference (such as microwave, high-frequency equipment operating environment), which may interfere with the product operation, making it unable to achieve the desired therapeutic effect.
- This device has not been tested for use in a magnetic resonance (MR) environment and should not be exposed to MR environments while patients are wearing the device. Keep it outside the magnetic resonance imaging (MRI) scanner room.
- Use the device according to the doctor or user manual's instruction.

Technical Information

Model/type	DVT001	Weight	355±10g
Battery Power	Powered by an internal 3.7V rechargeable lithium-ion battery	Degree of protection against electric shock	Type BF applied part
Air Pressure	0 – 55 mmHg (error: +/- 5 mmHg)	Type of protection against electric shock	Class II and built-in battery
Warranty	30 days	Grade of Ingress Protection	IP22
Recommend Treatment time	Prescribed by a physician	Treatment position	Calf
Air Bladder Chamber	Single bladder with a divided chamber	Power Charger	Input: AC100-240V, Output: DC5.0V-2A
Modes	1 mode	Mode of operation	Continuous operation
Software version	A0	Program Name	Auto Mode
The time required for the equipment to warm from the minimum storage temperature between uses until it is ready for intended use	30 minutes	The time required for the equipment to cool from the maximum storage temperature between uses until it is ready for intended use	15 minutes
Wrap cuff size (±10mm)	L: 600mm; W: 360mm	Main unit size (±1mm)	136X73X30mm
Typical service life of Battery	300 recharge cycles	Typical operation time after fully charged	Approx. 20 hours
Typical charging time	Approx. 4 hours		
Note: Device and its accessories are not intended to be sterilized.			
Not for use in an oxygen rich environment (>22% O ₂).			


Condition of Working, Storage and Transportation

- Normal working ambient temperature: 41 - 104°F (5 - 40°C)
- Normal working ambient humidity: 15 - 90% RH
- Store and transport ambient temperature: -13 - 158°F (-25 - 70°C)
- Store and transport ambient humidity: 0 - 90% RH
- Atmospheric pressure: 70 - 106kPa

Wrap Application

Wrap the leg wrap cuff around your calf and secure the strap to hold it in place. Ensure the wrap is secure, but not overly tight.


Note: The wrap may be placed directly against the skin or over a light dressing.

Power On

Hold the On/Off button for approximately 1.5 seconds to turn on the device. The buzzer will sound once when the device is powered on.

Upon start-up, the time display shows the accumulated working time, the air pressure display area shows 00, and the green battery indicators light up according to the battery level.


Working Mode

Once turned on, the device first displays the accumulated working time for 2 seconds. Afterward, the device starts inflating, and the display shows the real-time inflation pressure. Upon reaching the maximum pressure of 55mmHg, the device remains at 55mmHg for 4 seconds and then starts deflating. During deflation, the display of pressure is changed to "00", and keeps flashing until the next cycle starts. (Note: A complete circle of inflation, deflation, and rest takes approximately 60 seconds).

Power Off

Hold the On/Off button for approximately 1.5 seconds to turn off the device. The buzzer will sound twice when the device is shut down, the digital display turns off, and the green battery indicators turn off.

Charging

Use the provided power charger. The green battery indicators will keep flashing when the device is being charged, and become solid when the device is fully charged.

Note: The device can work normally during the period of charging. Use only the charger provided by manufacturer. Use of the wrong charger can cause damage and even danger to the device and patient.

The four green lights are used to represent the battery level of the device: Four green lights mean that the device is fully charged; Three lights mean that it has sufficient power; Two lights mean that the battery is consumed to the medium level; The battery level is low when there is one light.


When the battery power drops below a critical level, the low battery alert will occur. See section for Alerts.

Alerts

Low Battery:


When the battery is close to depleted, the device stops working, the first battery light keeps flashing quickly, the time display area shows "E2", the air pressure display area shows "LB", and the buzzer keeps beeping for 10 seconds. Afterward, the device turns off automatically. The device must be charged before next use.

Low pressure / Air leakage:

When the inflation time exceeds approximately 30 seconds and the inflated air pressure has not reached 55mmHg, the device stops working, all the battery lights keep flashing quickly, the time display area displays "E1", the pressure display area displays "LP", and the buzzer keeps beeping for 10 seconds. Afterward, the device automatically shuts down.

The wrap must be secured snugly around the calf to avoid low pressure alert. The device may be damaged if the wrap has been snugly applied and the following alerts appear: battery lights keep flashing, the display shows "LP", and the device continually shuts down. Please refer to the troubleshooting section.

Time Readout Display

Shows the cumulative number of hours the device has been in use. This number will remain on display when the device is on unless it is manually reset or if the hour count exceeds ninety-nine hours.


To reset the number manually:

- While the device is off, press and hold the power button for 10 seconds.
- The device will turn on and give an extended powering on sound.
- The total use time will be cleared and read "00".

Extender Cuff Application

Section A of the DVT unit may be cut for custom user fitting of the unit.

DO NOT CUT section B or any sewn seam as this may result in damage to the DVT unit.

Cleaning and Disinfecting

NOTE: Inspect the device and follow the cleaning and disinfecting procedures prior to each use. Cleaning and Disinfecting are recommended when required.

⚠ WARNING: Device must be turned off and disconnected from the wall outlet prior to and during cleaning or disinfecting. DO NOT IMMERSE DEVICE IN ANY LIQUID FOR ANY REASON.

- Clean the outer surface of the device using a soft cloth, moistened with soapy water or 70% isopropyl alcohol.
- Hand wash the exterior of the wraps using a soft cloth, moistened with soapy water or 70% isopropyl alcohol and let air dry.
- Do not use abrasive or volatile cleaners.
- To ensure the device is completely dry prior to use, leave the device in the OFF condition and disconnected from the wall outlet for at least 30 minutes after cleaning or disinfecting.
- Do not place device and the wrap in an autoclave, dryer or microwave.
- Do not use a hair dryer to accelerate drying.
- Do not place the device on top of or in front of portable or stationary radiators to accelerate drying.

Troubleshooting

There are no user serviceable parts. If your device is not operating properly, please check below for the common problems and suggested solutions. If the recommended action does not solve the problem, please contact Breg Customer Care at 1-800-321-0607 for assistance.

Problem	Possible Cause	Possible Solution
The skin turns red or the skin feels irritated.	The use time may be too long.	Reduce the use time.
No information showing on the device display.	The battery is depleted.	Recharge the battery.
Device stops operating during use.	The battery is depleted.	Recharge the battery.
Low pressure/air leakage alert triggers.	Leg wrap cuff may be applied too loosely.	Ensure the leg wrap cuff is applied snugly around the calf.

Maintenance and Disposal

- Inspect the device and all components for any damage that may have occurred prior to each use (for example, frayed or cut charging cord, cracked plastic housings, torn wraps, etc).
- Do not attempt to connect the power charger to a wall power supply if any damage is noticed.
- Avoid subjecting the device to shocks, such as dropping the device.
- Do not handle the leg wrap cuff with any sharp objects, to avoid damage to air bladder. If the air bladder is punctured or you notice a leak, do not attempt to repair the device or wrap.
- Avoid folding or creasing the air bladder during use and transportation of the device.
- Battery is not replaceable.
- If not in use for a long time, charge the battery every three months to prolong the battery life.

⚠ This device is an electromechanical device that includes printed circuit boards and rechargeable batteries. Dispose of the battery-containing device according to the local, state, or corresponding laws. Do not discard in landfill. Do not discard the device in regular waste.

Warranty

Breg, Inc. provides a limited warranty against manufacturing or materials defects for 30 days from initial purchase under normal use for which it was intended. The warranty applies to the device only, and the accessories are not covered by this warranty.

Breg's obligation under this warranty is limited to the replacement or repair of any defective part or parts of this product. No other warranty, expressed or implied, is given and no affirmation of or by seller, by words or action, will constitute a warranty.

For more information, please contact Breg Customer Care at 1-800-321-0607.

Electromagnetic Compatibility

⚠ WARNING: Don't use near active HF surgical equipment and the RF shielded room of an ME system for magnetic resonance imaging, where the intensity of EM disturbances is high.

⚠ WARNING: Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

⚠ WARNING: Use of accessories, transducers, and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

⚠ WARNING: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 12 in (30 cm) to any part of the equipment, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

All necessary instructions for maintaining BASIC SAFETY and ESSENTIAL PERFORMANCE regarding electromagnetic disturbances for the warranty period:

Guidance and manufacturer's declaration – electromagnetic emission			
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.			
Emission test	Compliance	Electromagnetic environment – guidance	
RF emissions CISPR 11.	Group 1	The device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.	
RF emissions CISPR 11.	Class B	The device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.	
Harmonic emissions IEC 61000-3-2	Class A		
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies.		
Guidance and manufacturer's declaration – electromagnetic immunity			
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±15 kV air	±8 kV contact ±15 kV air	Floors should be wood, concrete or ceramic tile. If floor is covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±2kV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	±1 kV differential mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% U_T (>95% dip in U_T) for 0.5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 sec	<5% U_T (>95% dip in U_T) for 0.5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the customer or the user of the device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery.
Power frequency (50Hz/60Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE: U_T is the a.c. mains voltage prior to application of the test level.			

Guidance and manufacturer's declaration – electromagnetic immunity

The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms	Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = 1,2\sqrt{P}$ $d = 1,2\sqrt{P}$ 80 MHz to 800MHz $d = 2,3\sqrt{P}$ 800 MHz to 2,5 Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ^a should be less than the compliance level in each frequency range. ^b
Radiated RF IEC 61000-4-3	10 V/m 80 MHz to 2.5 GHz	10 V/m	

NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a: Field strength from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the device.

b: Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Recommended separation distances between portable and mobile RF communications equipment and the device.

The device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)		
	150 kHz to 80 MHz	80 MHz to 800 MHz	800 MHz to 2.5 GHz
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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