
Polar Care® Kodiak®

Instructions for Use

Breg, Inc.
2382 Faraday Avenue, Suite 300
Carlsbad, CA 92008 U.S.A.
P: 800-321-0607
760-795-5440
F: 800-329-2734
www.breg.com
©2025 Breg, Inc.

AW-1.00564 Rev E 05/25

E/U authorized
representative
MDSS GmbH
Schiffgraben 41
D-30175 Hannover
Germany

UK Responsible Person
Joint Operations Ltd
11 Interface Business Park
Royal Wootton Bassett
Wiltshire
England
SN4 8SY

Breg Polar Care Kodiak

CAUTION: United States Federal Law restricts this device to sale by or on the order of a medical professional.

For questions or product inquiries contact **Breg customer care** at **1-800-231-0607** or **+1-760-795-5440**.

Table of Contents

Indications for Use	3
Symbol Key	3
Contraindications / Risk Factors	4
Warnings	6
Setup and Operation	7
Power Adapter	7
Materials	7
Troubleshooting	8
Electrical Safety	9
Emission and Immunity Table	10
Physical and Electrical Specifications	12
Warranty	12
Cold Therapy Protocol	13
Sample Cold Therapy Protocols	14
Patient Discharge Information Healthcare Professionals	16

Indications for Use

The extended use of motorized cold therapy modalities, such as Breg’s Polar Care Cold Therapy, have been shown to have many beneficial effects following surgery. Motorized cold therapy has repeatedly been shown to decrease postoperative pain, swelling, inflammation, and narcotic use following a variety of surgical procedures, such as surgery to the shoulder, knee, and back. For a bibliography of extended use motorized cold therapy studies, contact Breg, Inc. at 800-321-0607 or 760-795-5440.

Symbol Key

Symbols used on the Polar Care Unit

Symbol	Description
	This is the safety alert symbol. It is used to alert you to potential personal injury hazards. Obey all safety messages that follow this symbol to avoid possible injury or death.
	Type BF, applied part
	Class II Equipment
	Operating Instructions
	Refer to instruction manual / booklet
	Single patient multiple use.
IP22	Protected against solid foreign objects of 50 mm Ø and greater. Protected against vertically falling water drops when tilted up to 15°.
	Disposal of product according to local regulations.
	Medical Device

Contraindications / Risk Factors

Guidelines for Use

1. Screen Patients for Contraindications and Risk Factors

Before prescribing cold therapy, always consider the patient’s medical history, particularly any contraindications or risk factors. If not appropriately prescribed, motorized cold therapy can result in serious cold-induced injury, including full thickness skin necrosis.

CONTRAINDICATIONS	
Patients with any contraindications should not use Polar Care	
• History of cold injury, frostbite, or adverse reactions to local cold application. • Patients that are incoherent due to general anesthesia, sedation, or coma. • Application areas with compromised local circulation or potential wound healing problems, including localized compromise due to multiple surgical procedures.	• Circulatory syndromes, including Raynaud’s disease, Buerger’s disease, peripheral vascular disease, vasospastic disorders, sickle cell anemia, and hypercoagulable clotting disorders. • Local tissue infection. • Hand/wrist or feet/ankle surgery with polyneuropathy. • Diabetic Polyneuropathy.

RISK FACTORS for Cold-Induced Injury	
• Pathologic sensitivity to cold. • Behaviors that negatively effect circulation, including poor nutritional status, smoking and tobacco use, excessive caffeine use, and excessive alcohol use. • Patients with cold application area desensitization due to local anesthesia or regional nerve blocks. • Medications that have a negative effect on peripheral vascular circulation, including beta adrenergic blockers and local epinephrine use (such as in local anesthetics).	• Medications that have a negative effect upon mental capacity. • Excessive moisture at the application site due to excessive bleeding, sweating, or condensation. • Diabetes. • Hand/wrist or feet/ankle surgery. • Cognitive disabilities. • Communication barriers. • Young children and the elderly.

If the risk of cold-induced injury outweighs the benefits of cold therapy, do not prescribe Breg's Polar Care Cold Therapy. If you prescribe this product to patients with risk factors, consider taking special measures to control the risk, such as:

- Recommending more frequent skin checks.
- Requiring more frequent follow-up examinations.
- Using a thicker insulation barrier between the pad and skin.
- Prescribing shorter durations of application, less frequent application, or eliminating nighttime application.

2. Apply Insulation Barrier and Cold Therapy Pad

Always use an insulation barrier (such as Webril, Kerlix, cast padding, elastic bandage) between the Cold Therapy Pad and skin. **Do not let any part of the Pad touch skin.** If a sterile dressing has been applied to the treatment site that does not completely cover the skin under the pad, use an additional insulation barrier. The Pad alone can be too cold to be applied directly to the skin.

Use only Breg Cold Therapy Pads. Other pads may be colder, increasing the risk of skin injury.

Use only Breg Cold Therapy Pads designed for the body part being treated. Multi-Use and Rectangle Pads may be used on shoulders, knees, back, and hips. If the Cold Therapy Pad is placed on a body part that it is not designed for, or if another brand of pad is used, the skin can get too cold and be injured.

Do not cover the Cold Therapy Pad with dressing, wrapping, bracing or casting that prevents the patient from checking the skin under the Pad.

3. Provide Prescription and Instructions

Breg Polar Care Cold Therapy units are classified by the FDA as Class II medical devices that must be prescribed by a physician or licensed healthcare practitioner. A proper prescription for use must include:

- Frequency and duration of use (and breaks if applicable);
- Frequency and instruction on skin inspections;
- Treatment Period.

Use the Physician Prescription Form on the other side of this document. This document should be given to the patient (or caregiver) upon discharge or transfer from the recovery room. Common operating temperatures for patients without risk factors for extended use (longer than 20 minutes) are 45° to 60° F (7° to 16° C) for core regions and 50° to 60° F (10° to 16° C) for extremities. When used with a body part appropriate Cold Therapy Pad and an insulation barrier, the average operating temperature of the Breg Polar Care Cold Therapy System will automatically fall within these ranges.

A variety of extended use motorized cold therapy protocols are reported in the literature. Treatment protocols vary based upon specific patient conditions and health history, physician experience with motorized cold therapy application, and medical judgment. See the Sample Cold Therapy Protocols to the right.

Tell patients how to inspect skin under the Cold Therapy Pad without compromising the sterile site and how often to do so.

Instruct patients to stop using Breg's Polar Care Cold Therapy and contact you immediately if they experience any adverse reactions such as: increased pain, burning, increased swelling, itching, blisters, increased redness, discoloration, welts, or other changes in skin appearance.

Review all product information with your patient.

Review the Information for Patients and Medical Professionals in this document, the Operating Instructions on the reverse side of this product insert, and the Cold Therapy Pad Fitting Instructions.

Warnings

Information for Patients & Medical Professionals

WARNINGS

Breg's Polar Care Cold Therapy can be cold enough to seriously injure skin. Follow this information, the Operating Instructions in this product insert, and the Cold Therapy Pad Fitting Instructions (provided with each Cold Therapy Pad).

Discuss Treatment with your Licensed Health Care Practitioner


Provide a complete medical history including any reactions to cold. Certain medical conditions make cold-induced injury more likely. Ask your practitioner about potential adverse reactions and cold-induced injuries.

Use Only as Prescribed

Use only as instructed by your practitioner regarding the frequency and duration of cold application and breaks, how and when to inspect the skin, and total length of treatment.

Skin Checks

Inspect the skin under the pad (by lifting the edge) as prescribed, typically every 1 to 2 hours. Do not use if dressing, wrapping, bracing, or casting over the pad prevents skin checks. Stop using and contact your practitioner immediately if you experience any adverse reactions, such as increased pain, burning, increased swelling, itching, blisters, increased redness, discoloration, welts, other changes in skin appearance, or any other reaction identified by your practitioner.

General Safety


Keep cords and hoses away from the neck to avoid risk of strangulation. No modification of this equipment is allowed.

The device is intended for single patient multiple use.

Apply Insulation Barrier & Intelli-Flo® Pad

Always use an insulation barrier (such as Webril, Kerlix, cast padding or elastic bandage) between the Intelli-Flo Pad and skin. Do not let any part of the pad touch skin. If a sterile dressing has been applied to the treatment site that does not completely cover the skin under the pad, use an additional insulation barrier. Use only Breg Intelli-Flo Pads designed for the specific body part. Other pads may be colder, increasing the risk of skin injury.

1. With the insulation barrier in place, apply the pad with the blue side up.
2. Use elastic straps to hold the pad in place.


Notice to Users

If any serious incident has occurred in relation to the device, report to Breg and the competent authority of the Member State in which the user and/or patient is established.

Setup and Operation


1. To fill:

- Unlock handle and remove lid.
- Fill with cold water to line, then ice to line.
- Replace lid and lock with handle in upright position.

2. To connect unit to pad:

- Push down on the two metal tabs.
- Push connectors together until metal tabs pop up.


3. To turn on/off:

- Plug power supply into wall and receptacle on lid.
- Operation begins automatically with power.
- Unplug power supply to turn off.


4. To refill:

- Turn off.
- Press down on two metal tabs to disconnect unit and pad.
- Remove lid and place upside down on a flat surface.
- Drain and follow fill instructions.

Power Adapter

- Align the lock to the release tab
- Press the plug and turn clockwise
- To unplug, press the release button and rotate the plug in the counterclockwise direction

Materials

Polar Care Kodiak:

- PVC
- Nylon
- EPDM
- Silicone
- Nitrile
- Polypropylene
- Polyether
- Copper

Intelli-Flo Pads:

- Tricot/Ply Blend Foam
- Dureflex Poly Film
- EPDM rubber
- Polypropylene

Troubleshooting

Problem: Pump not running, water not flowing to pad, or pad not cold

Possible solutions:

- Use larger ice for optimal performance.
- Allow 10 minutes for flow and pressure to stabilize.
- Ensure power outlet is working and plugs are fully engaged.
- Ensure unit has both ice and water.
- Ensure hose is not kinked.
- Disconnect and reconnect the pad and unit.
- Release air by depressing the black plus-shaped part inside the unit connector. Note: water may be released.
- Place unit on a table or other raised surface.
- Decrease tension of bandages or straps around the pad.
- Remove pad and lay it flat. Allow pad to fill; reapply.
- Clean filter: Disconnect pad. Remove unit lid. Pull filter cap from bottom of lid. Remove foam filter. Rinse filter cap and filter to remove clogs. Reassemble.

Problem: Condensation

Possible solutions:

- Wrap material over pad and hose to minimize air exposure.
- Protect the wound site by using a sterile dressing with waterproof barrier.

Problem: Unit is leaking

Possible solutions:

- Disconnect unit connector. Ensure the metal tabs are pressed down; reconnect.
- Apply lubricant to o-rings of connectors.
- If leaking continues, or if a leak is detected in the pad or unit lid, stop using the unit and contact Breg Customer Care at 800-321-0607 or 760-795-5440.

Cleaning, Maintenance & Service

After use, drain and dry pump with a soft cloth. Warm water and mild detergent may be used occasionally to clean the pump and tubes. There are no serviceable parts. Contact Breg Customer Care for replacement parts.

Disposal

Dispose of product according to local regulations.

Electrical Safety

Like all electrical products, you must follow precautions to avoid electrocution, fire, burns, or other injuries. Operator do not touch exposed transformer pins and patient at the same time. Use only the transformer provided with the Breg Polar Care Cold Therapy unit to meet the requirements of EN 60601-1.

Keep electrical connections dry and free of debris

Even though the Polar Care unit is designed to hold and pump water, the electrical connections at the power outlet and from the power cord to the Polar Care unit must be kept dry. Do not handle the transformer or electrical cord with wet hands. Always keep the unit in a place where the connections will not fall into water (e.g., a tub, sink, etc.). If the electrical connections fall into water, do not touch any wet part of the product. Unplug only at dry electrical connections.

Avoid flammables and oxidizers

Do not use in places with flammable vapors or gasses (e.g. flammable anesthetics), high oxygen concentrations, or other oxidizers (e.g. nitrous oxide).

Protect the power cord and transformer

Keep the cord and transformer away from heated surfaces. Set up and use this equipment in a low traffic location away from children and pets. Never operate this product if it has a damaged cord, transformer, or plug, which could result in electric shock. If the cord, transformer, or plug is damaged, unplug and contact Breg Customer Care, (800) 321-0607 or (760) 795-5440.

Unplug to stop the product when not in use

Unplugging the Polar Care unit turns it off. Always unplug immediately after use. Never leave it plugged in while unattended. Never run pump without water. Running the unit without water will cause permanent damage.

Electromagnetic Interference and Emissions

This Polar Care unit may cause electromagnetic or other interference with other electrical devices. To check whether the Polar Care unit is interfering with another device, unplug the Polar Care unit. If this corrects the problem, move the Polar Care unit or other device, or use an outlet on a different circuit.

Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation.

Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the Kodiak including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

Emission and Immunity Table

Electromagnetic Emissions Declaration Table		
The Kodiak is suitable for use in the specified electromagnetic environment. The consumer and/or user of the Kodiak should ensure that it is used in an electromagnetic environment as described below:		
Emission Test RF Emissions CISPR 11	Compliance Group 1	Electromagnetic environment - guidance The Kodiak uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause interference in nearby electronic equipment.
RF Emissions CISPR 11	Class B	The Kodiak is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic Emissions IEC 61000-3-2	Class B	
Voltage fluctuations / flicker emissions IEC 61000-3-3	Complies	

Immunity Test	IEC 60601-1-2 Test Level	Electromagnetic Environment Guidance
Conducted Disturbances induced by RF Fields IEC 61000-4-6	3V 0.15MHz – 80MHz 6V in ISM & amateur radio bands between 0.15Mz – 80MHz 80% AM at 1kHz	Portable and mobile RF communications equipment should be used no closer to any part of the Kodiak, including cables, than the recommended separation distance calculated from the equation appropriate for the frequency of the transmitter. Recommended Separation Distance $d=1,2\sqrt{P}$ $d=1,2\sqrt{P}$ 80 MHz to 800 MHz $d=2,3\sqrt{P}$ 800MHz to 2,5 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey (a), should be less than the compliance level in each frequency range (b). Interference may occur in the vicinity of equipment marked with the following symbol:
Radiated RF EM fields IEC 61000-4-3	10 V/m 80 MHz – 2.7 GHz 80% AM at 1 kHz	

Electromagnetic Immunity Declaration Table

The Kodiak is suitable for use in the specified electromagnetic environment. The user of the Kodiak should ensure that it is used in an electromagnetic environment as described below:

Immunity Test	IEC 60601-1-2 Test Level	Electromagnetic Environment Guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV 100 kHz repetition frequency	Main power quality should be that of a typical commercial and/or hospital environment.
Surges (line to line) IEC 61000-4-5	± 0,5 kV, ± 1 kV	Main power quality should be that of a typical commercial and/or hospital environment.
Proximity Fields from RF Wireless Communications Equipment IEC 61000-4-3	Reference IEC 60601-1-2:2014+AMD1:2020 table 9	
Voltage Dips and Fluctuations IEC 61000-4-11	0% for ½-cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270°, and 315° 0% for 1-cycle and 70% for 30 cycles at 0°	Main power quality should be that of a typical commercial and/or hospital environment. If the user of the Kodiak requires continued operation during power main interruptions, it is recommended that the Kodiak be powered from an uninterruptible power supply or a battery.
Voltage Interruptions IEC 61000-4-11	0% for 300 cycles	
Rated Power Frequency Magnetic Fields IEC 61000-4-8	30 A/m 50Hz	
		If improper operation occurs, it may be necessary to position the Kodiak further from sources of power frequency magnetic fields or to install magnetic shielding. The power frequency magnetic field should be measured in the intended installation location to ensure that it is sufficiently low.

Proximity Magnetic Fields IEC 61000-4-39	IMMUNITY TEST LEVELS TO DISTURBANCES FROM PROXIMITY MAGNETIC FIELDS		
	Test Frequency	Modulation	Immunity Test Level (A/m)
	30 kHz	CW	8
	134.2 kHz	Pulse modulation ¹ 2.1 kHz	65
	13.56 MHz	Pulse modulation ¹	7.5
		50kHz	

¹Modulate carrier using 50% duty cycle square wave. Immunity test levels specified as RMS before modulation is applied.

Physical and Electrical Specifications

Unit Rated: 6Vdc, 600-1000mA Power supply is specified as part of ME Equipment	
Manufacturer: MERRYKING ENTERPRISES (HK) CO., LTD. or Group West Rated Input: 100-240V~, 50-60Hz	Model: MKB2-0601000DEXD or 18UR-06-1000 Output: 6Vdc, max. 1.0A
Environmental Requirements: Transport & Storage: Operating:	-25°C (-13°F) no humidity control to 70°C (158°F) up to 93%RH 5°C to 40°C (41°F to 100°F)15-93%RH 700 hPa to 1060 hPa
Standards Compliance:	Designed to conform to applicable requirements of: EN 60601-1, 3rd Edition
Electrical Classifications:	Class II, Continuous Operation. Not suitable with flammable anesthetics.

Product #:	Weight	Average Operating Range	Expected Service Life
100291-000	4.0 lbs (1.8kg) Empty 14.5 lbs (6.6kg) Filled	Average operating range of Intelli-Flo Pad with Polar Dressing ≥ 45°F (≥7.2°C).	1 year from date of first use.

Warranty

Breg, Inc. warrants that this product is free from defects in workmanship, materials, and fitness for use for 180 days from initial purchase under normal use for which it was intended and under direct supervision of a licensed health care practitioner. Breg, Inc.'s obligation under this warranty is limited to the replacement or repair of any defective part or parts of this product. All expressed or implied warranties, including the warranty of merchantability and fitness for a particular purpose, are limited to the actual warranty period set forth above. No other warranty, expressed or implied, is given and no affirmation of or by seller, by words or action, will constitute a warranty.

Cold Therapy Protocol*

Treatment Period	Awake / Asleep	Frequency / Duration	Inspect Skin Every
Day:	Awake		
Through			
Day:	Asleep		
Day:	Awake		
Through			
Day:	Asleep		
Day:	Awake		
Through			
Day:	Asleep		

*To be completed by a licensed healthcare professional.

Discuss Treatment with your Licensed Health Care Practitioner

Provide a complete medical history including any reactions to cold. Certain medical conditions make cold-induced injury more likely. Ask your practitioner about potential adverse reactions and cold-induced injuries.

Use Only as Prescribed

Use only as instructed by your practitioner regarding the frequency and duration of cold application and breaks, how and when to inspect the skin, and total length of treatment.

Skin Checks

Inspect the skin under the pad (by lifting the edge) as prescribed, typically every 1 to 2 hours. Do not use if dressing, wrapping, bracing, or casting over the pad prevents skin checks. Stop using and contact your practitioner immediately if you experience any adverse reactions, such as increased pain, burning, increased swelling, itching, blisters, increased redness, discoloration, welts, other changes in skin appearance, or any other reaction identified by your practitioner.

Sample Cold Therapy Protocols

The following protocols are examples of appropriate post-operative extended (>20 minute application) motorized cold therapy for patients with no risk factors. Therefore, these sample protocols would not apply to use on extremities, which may require extra precaution as described in Section 1. This is not a comprehensive list of treatment possibilities. Treatment protocols for extended motorized cold therapy may vary based upon specific patient conditions and health history, physician experience with cold therapy, and medical judgment.

When used with a body part appropriate Cold Therapy Pad and an insulation barrier, the average operating temperature of the Breg Polar Care Cold Therapy System will automatically fall within these ranges.

Breg Polar Care Cold Therapy products have the ability to deliver cold after ice replacement for approximately 6 to 8 hours.

The duration of a cyclic application may vary depending upon the patient. If the patient does not experience pain relief, the physician may increase the duration of application. As the application duration is increased, the frequency of the skin inspections should increase.

Sample Protocol 1

Treatment Period	Frequency /Duration	Temp.	Skin Every
Day 1-3	While awake: Continuous	45° - 60° F (7° - 16° C)	Inspect skin under pad every 1-2 hours
	While asleep: Continuous	45° - 60° F (7° - 16° C)	Upon waking
Day 4-10	While awake: Cyclic: 1 hour on and 1 hour off	45° - 60° F (7° - 16° C)	Periodically inspect the skin under pad
	While asleep: Continuous	45° - 60° F (7° - 16° C)	Upon waking
Day 11 and beyond	While awake: As needed for pain control: Continuous for 1 hour intervals; not to exceed 12 hours/day	45° - 60° F (7° - 16° C)	Periodically inspect the skin under pad

*50° - 60° F (10° - 16° C) for pads intended for use on hands, feet or face.

Sample Protocol 2

Treatment Period	Frequency /Duration	Temp.	Skin Every
Day 1-3	While awake: Continuous	45° - 60° F (7° - 16° C)	Inspect skin under pad every 1-2 hours
Day 4 to suture removal (Day 7-14)	While awake: As needed for pain control: Continuous for 1 hour intervals; not to exceed 12 hours/day	45° - 60° F (7° - 16° C)	Periodically inspect the skin under pad

Sample Protocol 3

Treatment Period	Frequency /Duration	Temp.	Skin Every
Day 1-3	While awake: Cyclic: 1 hour on and 1 hour o	45° - 60° F (7° - 16° C)	Periodically inspect the skin under pad
	While asleep: Continuous	45° - 60° F (7° - 16° C)	Upon waking
Day 4 and beyond	As needed for pain control: Continuous for 1 hour intervals; not to exceed 12 hours/day	45° - 60° F (7° - 16° C)	Periodically inspect the skin under pad

Patient Discharge Information Healthcare Professionals

Follow this protocol prior to discharging the patient from facility care to home use.

1. **Patient Screen.** Screen the patient for any contraindications and/or associated risk factors. If the patient has any contraindications, do not dispense Breg's Polar Care Cold Therapy to the patient. If the patient has any associated risk factors, consult with the licensed healthcare practitioner to determine the appropriateness of application of Breg's Polar Care Cold Therapy to that patient.
2. **Instructions For Use.** Instruct the patient on how to properly use Breg's Polar Care Cold Therapy. Review the Operating Instructions affixed to the unit with each patient.
3. **Prescription.** Instruct the patient regarding the licensed healthcare practitioner's prescribed protocol (frequency and duration of use and breaks if applicable), frequency and instruction on skin inspections, and treatment period.
4. **Potential For Injury.** Inform the patient that improper use can result in serious skin injury, including necrosis. Emphasize the importance of following the prescribed protocol, proper pad application, and skin inspection.
5. **Proper Pad Application.** Instruct the patient that an insulation barrier must be between the Cold Therapy Pad and skin during use. No part of the Cold Therapy Pad should touch the skin for any period of time.
6. **Skin Inspection.** Instruct the patient to inspect the skin receiving cold treatment per the practitioner's instructions, typically every 1 to 2 hours. If dressing, wrapping, bracing, or casting over the Cold Therapy Pad prevents the ability of the patient to regularly perform skin checks under the pad, do not dispense Breg's Polar Care Cold Therapy to the patient.
7. **Discontinue.** Instruct the patient to stop using Breg's Polar Care Cold Therapy and contact their licensed health care practitioner immediately if they experience any adverse reactions such as: increased pain, burning, increased swelling, itching, blisters, increased redness, discoloration, welts, or other changes in skin appearance.
8. **Documentation.** Give the patient this document with the Cold Therapy Protocol (on reverse side, top left) filled out.

Breg Polar Care Cold Therapy products have the ability to deliver cold after ice replacement for approximately 6 to 8 hours.

The duration of a cyclic application may vary depending upon the patient. If the patient does not experience pain relief, the physician may increase the duration of application. As the application duration is increased, the frequency of the skin inspections should increase.

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