ENGLISH

INTENDED USE/INDICATIONS:

Indicated for dislocations, clavicle and humerus fractures.

INTENDED USERS:

The intended user should be a licensed professional, the patient, or patient’s caregiver. The user should be able to read, understand and be physically capable to perform all the directions, warnings and cautions provided in the information for use.

CONTRAINDICATIONS:

Not Applicable.

⚠️ WARNINGS AND PRECAUTIONS:

- Carefully read fitting instructions and warnings prior to use. To ensure proper performance of the brace, follow all instructions.
- If you experience, or have an increase in pain, swelling, skin irritation, or any adverse reactions while using this product, immediately consult your medical professional.
- This device will not prevent or reduce all injuries. Proper rehabilitation and activity modification are also an essential part of a safe treatment program. Consult your licensed health care professional regarding safe and appropriate activity level while wearing this device.
- Follow application and care process described in instructions.
- For single patient use only.

INSTRUCTIONS FOR USE:

- Position sling against rib cage. Bring web straps straight up and over the patient’s shoulders. Place affected arm through the sling.
- Criss cross both straps against the patient’s back.
- Attach hook straps to the front of the sling in a criss cross motion supporting the affected arm and adjust to length.
- Secure arm against the body by attaching the swathe strap to the sling. Place one end of the strap to the sling. Encircle the affected arm and bring around the patient’s back and under the non-affected arm. Secure loose end to the front of the sling.

USE AND CARE:

Fasten hook and loop straps. Hand wash in 30°C (86°F) water with mild soap, rinse thoroughly, Rinse and lay flat to air dry. Do not machine wash or tumble dry.

MATERIAL CONTENT:

Acetal, Cotton, Nylon, Polyester, Polyurethane

WARRANTY:

Breg products come with a limited warranty against manufacturing and material defects. No other warranty, expressed or implied, is given and no affirmation of or by seller, by words or action, will constitute a warranty. Please see www.breg.com

DISPOSAL:

Dispose of product according to local regulations.

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