ENGLISH

INTENDED USE/INDICATIONS:

Soft knee bracing products are designed to provide control and support for ligament injuries and instabilities. Indicated for most ACL, PCL, MCL, LCL injuries and instabilities, meniscal injuries, mild OA, Medial and Lateral instabilities.

INTENDED USERS:

The intended user should be a licensed professional, the patient, or patient’s caregiver. The user should be able to read, understand and be physically capable to perform all the directions, warnings and cautions provided in the information for use.

CONTRAINDICATIONS:

Not Applicable.

⚠️ WARNINGS AND PRECAUTIONS:

- This device is intended for use only under the direction of a medical professional. If you experience any pain, swelling, sensation changes or any unusual reactions while using this product, consult your medical professional immediately.
- This product contains natural rubber latex, which may cause allergic reactions.
- For single patient use only.

INSTRUCTIONS FOR USE:

- Pull up the sleeve over the knee joint until the patella is centered in the middle of the circle pattern.
- Tighten the lower strap first, then the upper strap on the back of the leg, followed by the upper strap on the front of the leg (Figure A & B).
- Align the axis on the hinges with the center of the knee joint. Make strap adjustments if necessary (Figure C).

USE AND CARE:

Hand wash in 30°C (86°F) water with mild soap, rinse thoroughly. Rinse and lay flat to air dry. Do not machine wash or tumble dry.

MATERIAL CONTENT:

Natural Rubber Latex, Nylon, Polyester, Polyurethane, Spandex

WARRANTY:

Breg products come with a limited warranty against manufacturing and material defects. No other warranty, expressed or implied, is given and no affirmation of or by seller, by words or action, will constitute a warranty. Please see www.breg.com

DISPOSAL:

Dispose of product according to local regulations.

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