

Genesis Walker Boot



The lightest full shell walker boot on the market

Walker / Ankle / Foot 



Just
Got
Better!



Genesis Walker Tall & Mid-Calf


Genesis Walker

Not all walker boots are created equal.

The Genesis Walker, inspired by **athletic footwear technology**, offers a **biomechanically natural gait**, mirroring a normal shoe for **seamless daily life integration** during recovery.



The lightweight, low profile **Genesis Walker** revolutionizes patient **comfort, maneuverability and control**.



The Genesis Advantage

Comfort

Comfort with Every Step

The **Genesis Walker is the lightest full shell pneumatic walker on the market***

It has been computer engineered to provide strength where required while eliminating extra weight.



Dynamic Air Stride Technology
Inspired by athletic shoes, the Genesis Walker features a **heel air chamber** to **absorb shock** and enhance walking comfort.

Customizable Shell

The Genesis Walker may be adjusted for patients with larger calf muscles. A cut line allows the boot to be cut to accommodate a large calf.




Happy feet, Healthy body

Narrow Heel Design
Mimics standard shoe width to reduce knee strain and facilitate easier, more **natural walking**.



Proper heel height and width reduces pelvic tilt and secondary knee pain

Low Heel Advantage
The Genesis Walker features the market's **lowest heel height**, aligning closely with standard shoes to minimize pelvic tilt and **prevent secondary hip and back pain**.



Maneuverability

* Data on file.

** Standard shoe is defined as a standard men's dress shoe with an approximate heel height of 1 inch.

The Genesis Advantage


Walk Normally While Recovering

The Genesis Walker's **streamlined design** provides a low profile look and feel. This reduces the feeling of bulk and makes it noticeably easier to maneuver.

Engineered Tread Design

Rocker Bottom

Engineered for natural gait by **mimicking healthy ankle movement**, promoting smoother walking with an immobilized ankle.



Scalloped No-Slip Tread

Athletic shoe tread design provides traction on a variety of surfaces permitting better cross-directional maneuverability.


Control

Performance that Delivers

The Genesis Walker delivers control, support and protection to the foot and ankle from time of injury until full recovery.


Dual Air Chambers

Offer **uniform compression** from foot to ankle, securing position for **enhanced control and comfort during healing**.



Added Cushioning in Air Bladder

Superior ankle support and comfort.



Built-in Pneumatic System

Integrated Air Pump

Enables **precise pressure adjustments** for foot and ankle support without the need for an external hand pump.

Hard Shell Convenience

Offers **cast-like control** with walker flexibility, adjustable for swelling, easy inspections, and dressing changes, ensuring optimal injury protection.



Enhanced Comfort

Removing binding from ankle and **Arser foot pad**.

Genesis Mid-Calf Walker

Utilizing the same great features as the Genesis Walker, the Genesis Mid-Calf Walker offers maximum **comfort**, **manueverability**, and **control** in a more compact size.

Great things come in **small packages**

Using the same cutting-edge design as the tall Genesis Walker, the mid-calf boot's short stature allows it to provide **optimal support** with **minimal impact** to your life.



Mid-Calf Design
At 25% shorter than the standard Genesis Walker, the Genesis Mid-Calf is low profile and ultra lightweight.

The lightest full shell walker boot just got lighter



Available with **custom logos**

- Genesis Walker Tall
- Genesis Walker Mid-Calf

Genesis Updated Features & Benefits



Introducing the enhanced Genesis Walker, meticulously designed for ultimate comfort, customization, and lightweight mobility. Revel in unparalleled plush comfort with our new luxurious foam cushioning and soft, trimmable toe cover for a snug fit. Benefit from increased durability and strength, thanks to our reinforced straps and superior stitching. Our extended forefoot pad ensures comprehensive comfort from above the ankle to the toes. Embrace the Genesis Walker's upgrades for superior, lightweight support and an elevated experience.

- **Added Cushioning:** Luxurious foam cushioning for superior ankle support and plush comfort.

- **Enhanced Comfort:** Removal of the binding from the ankle area eliminates unnecessary pressure points on front of foot.

- **Soft toe cover:** The soft toe cover is removable and trimmable.



- **Stronger Straps:** Reinforced stitching and overlapped strap material for superior strength.

- **Larger Foot Pad:** The forefoot pad now extends from above the ankle to the toes, dramatically enhancing patient comfort.

Genesis Part Numbers & Pricing

Description	Part Numbers	HCPCS	Logo
Mid-Calf			
WALKER GENESIS MID-CALF BOOT XS	BL525001	L4361OTS or L4360CF	Yes
WALKER GENESIS MID-CALF BOOT S	BL525003	L4361OTS or L4360CF	Yes
WALKER GENESIS MID-CALF BOOT L	BL525007	L4361OTS or L4360CF	Yes
WALKER GENESIS MID-CALF BOOT XL	BL525009	L4361OTS or L4360CF	Yes
WALKER GENESIS MID-CALF BOOT M	BL525005	L4361OTS or L4360CF	Yes
3 Straps			
WALKER GENESIS 3 STRP GRAY XS	BL515001	L4361OTS or L4360CF	Yes
WALKER GENESIS 3 STRP GRAY S	BL515003	L4361OTS or L4360CF	Yes
WALKER GENESIS 3 STRP GRAY M	BL515005	L4361OTS or L4360CF	Yes
WALKER GENESIS 3 STRP GRAY L	BL515007	L4361OTS or L4360CF	Yes
WALKER GENESIS 3 STRP GRAY XL	BL515001	L4361OTS or L4360CF	Yes
4 Straps			
WALKER GENESIS 4 STRP GRAY XS	BL514001	L4361OTS or L4360CF	Yes
WALKER GENESIS 4 STRP GRAY S	BL514003	L4361OTS or L4360CF	Yes
WALKER GENESIS 4 STRP GRAY M	BL514005	L4361OTS or L4360CF	Yes
WALKER GENESIS 4 STRP GRAY L	BL514007	L4361OTS or L4360CF	Yes
WALKER GENESIS 4 STRP GRAY XL	BL514009	L4361OTS or L4360CF	Yes
Accessories			
KIT, GENESIS COOL LINER, M	100507-030	L4392	No
KIT, GENESIS COOL LINER, L	100507-040	L4392	No
KIT, GENESIS COOL LINER, XL	100507-050	L4392	No
KIT, GENESIS COOL LINER MC, XS	100508-010	L4392	No
KIT, GENESIS COOL LINER, XS	100507-010	L4392	No
KIT BOOTIE/STRAP GENESIS XS	BL515041	L4002 and L4392	No
KIT BOOTIE/STRAP GENESIS S	BL515043	L4002 and L4392	No
KIT BOOTIE/STRAP GENESIS M	BL515045	L4002 and L4392	No
KIT BOOTIE/STRAP GENESIS L	BL515047	L4002 and L4392	No
KIT BOOTIE/STRAP GENESIS XL	BL515049	L4002 and L4392	No
KIT STRAP GENESIS XS	BL515071	L4002	No
KIT STRAP GENESIS S	BL515073	L4002	No
KIT STRAP GENESIS M	BL515075	L4002	No
KIT STRAP GENESIS L	BL515077	L4002	No
KIT STRAP GENESIS XL	BL515079	L4002	No
EXTENSION, WRAP, GENESIS	100692-000	L9900	No
ACCSY KIT, ACHILLES WEDGE, GENESIS S-M	100566-010	L9900	No
ACCSY KIT, ACHILLES WEDGE, GENESIS L-XL	100566-030	L9900	No
KIT HARD TOE GUARD GEN XS	BL515061	L9900	No
KIT HARD TOE GUARD GEN S	BL515063	L9900	No
KIT HARD TOE GUARD GEN M	BL515065	L9900	No
KIT HARD TOE GUARD GEN L	BL515067	L9900	No
KIT HARD TOE GUARD GEN XL	BL515069	L9900	No

