



Pinnacle Lumbar Brace Family

Introducing Breg's premium lumbar braces

The Pinnacle family of spinal orthoses considers patient comfort in every detail. These braces offer superior support, clinically advantageous features such as incision openings and ergonomic handles, and they are enhanced with premium materials and finishes.

Pinnacle is the peak of comfort for spinal recovery.

All of our lumbar braces have patient-centric design features, including:

- Lightweight and breathable Airmesh[®] padding
- Universal sizing that makes them easy to fit
- 3:1 mechanical advantage lacer system
- Dual-pull lacer system with independent top and bottom adjustment
- Ergonomic pull handles

Pinnacle Lumbar Orthoses

Breg's new Pinnacle lumbar braces are designed for patients who want support throughout recovery and those who want a sleeker look as they return to daily activities.



Pinnacle LO 627/642

Ergonomic pull tab
(on all Pinnacle braces)



Pinnacle LSO 631/648

Additional support up to T9 vertebra



Features of our Pinnacle lumbar orthoses include:



Pinnacle LSO 637/650

Additional lateral supports



Pinnacle LSO 639/651

Adjustable height xiphoid extension




Common examples of use for lumbar orthoses

- Post-operative Support
- Fracture management
- Post Laminectomy Syndrome
- Spondylolisthesis
- Decompressive Procedures
- Chronic Low Back Pain
- Deconditioned Trunk Musculature
- Spinal Stenosis
- Sprains / Strains
- Radiculopathy
- Spondylosis
- Bulging / Herniated Disc
- Degenerative Disc Disease
- Severe DJD
- Facet Syndrome
- IDET Procedures
- SI Dysfunction
- Sacroiliitis


Pinnacle LO 627/642

SP20627-000



Pinnacle LSO 631/648

SP20631-000



+ Additional support up to T9 vertebra

Pinnacle LSO 637/650

SP20637-000



- + Additional lateral supports (provided with brace)

Pinnacle LSO 639/651

SP20639-000




- + Adjustable height xiphoid extension
- + Rigid anterior plate supports from symphysis pubis to xiphoid process



Pinnacle Thoraco-Lumbo-Sacral Orthoses (TLSO)

Breg's new Pinnacle TLSO brace is designed for patients who want support throughout recovery and those who want a sleeker look as they return to daily activities.



Features of our Pinnacle TLSO orthoses include:

Over-the-shoulder or under-arm strap configurations



Ability to convert to LO 627

High-end, simulated suede material for a polished aesthetic

No-tool-required anterior thoracic extension adjustment

Universal sizing



Easy-grip application strap

Pinnacle TLSO 464

SP20464-000



Common examples of use for TLSO braces

- Post-operative Support
- Fracture Management
- Decompressive Procedures
- Degenerative Disc Disease
- Kyphosis
- Scoliosis
- Kyphosis Secondary to Osteoporosis
- Spondylolisthesis
- Spinal Stenosis
- Spondylosis
- Bulging / Herniated Disc
- Severe DJD



Part Numbers

Lumbar Orthoses

Pinnacle LO 627/642	SP20627-000
Pinnacle LSO 631/648	SP20631-000
Pinnacle LSO 637/650	SP20637-000
Pinnacle LSO 639/651	SP20639-000
Pinnacle TLSO 464	SP20464-000

Accessories

Extension Panel, Pinnacle Belt	SP20010-000
Upgrade Kit, Pinnacle TLSO 464	SP21464-000
Pad Kit, Pinnacle Lumbar	SP20020-000



For more information visit
www.breg.com/lumbar

