



DVT Guardian™

Prevention | Portability | Dependability

Trusted DVT Prevention with Dual-Chamber Compression for Optimal Blood Flow

Introducing DVT Guardian your ideal Deep Vein Thrombosis (DVT) prevention partner for a safe and comfortable recovery at home or in a clinical setting. The innovative device uses dual-chamber compression to enhance blood circulation and aid in the prevention of DVT. With preset pressure cycles that simulate muscle contractions, the DVT Guardian offers up to 20 hours of operation on a single battery charge. Its portable, rechargeable design ensures you receive consistent protection.



Air Pressure, Time & Battery Level at a Glance



DUAL-CHAMBER
COMPRESSION



ONE-TOUCH
OPERATION



EASY-TO-READ
DISPLAY



CORDLESS
AND PORTABLE



TRIMMABLE
CUFF



UP TO 20 HOURS
BATTERY LIFE



LIGHTWEIGHT
WRAP



HOOK & LOOP
STRAPS



DVT Guardian™

Am I at Risk?

After surgery, patients can be at risk of developing deep vein thrombosis (DVT) due to reduced mobility and increased susceptibility to blood clotting. DVT occurs when a blood clot forms in a vein deep in the body, usually in the calf or thigh. If the clot dislodges, it can travel through the bloodstream and potentially block blood flow to the lungs, posing a serious health risk.*


DVT Guardian: Dependable DVT Prevention

Included Components: One pair of devices and a battery charger.

Device Display Features: Battery level indicator, air pressure display, and time display for easy monitoring.

Compression Cycle: One-minute cycle duration delivering 55mmHg of compression.

Cuff Extender: Optional Cuff Extender available to accommodate larger sizes.



□ INDICATIONS FOR USE:

The DVT Guardian can be used to:

- Aid in the prevention of Deep Vein Thrombosis
- Enhance blood circulation
- Diminish post-operative pain and swelling

PART NUMBERS

DVT001	DVT Guardian Device
DVT002	Charger, 2 Prong, 6 ft
DVT003	Charger, 2 Prong, 15 ft
DVT004	Cuff Extender, pair



Breg, Inc. 2382 Faraday Avenue, Ste. 300, Carlsbad, CA 92008 USA,
Phone: 800.321.0607, Fax: 760.795.5295, customer@BREG.com, BREG.com

©2025 BREG, Inc. All rights reserved. DVT Guardian is a trademark of BREG, Inc. For more information on Indications for Use, Contraindications, Warnings and Precautions please reference DVT Guardian Instructions for Use.

*Source: <https://orthoinfo.aaos.org/en/diseases--conditions/deep-vein-thrombosis/>