

Vectra<sup>®</sup> Pro

## A Legacy Reborn

Experience the evolution of  
reliability and innovation  
in recovery





# Introducing *Vectra Pro*

Experience the evolution of reliability and innovation in medical walking boots. Building on 17 years of trusted performance, the Vectra Pro raises the bar with advanced features for superior protection and comfort. Ideal for recovery from injuries or surgery, the Vectra Pro provides patients with exceptional support and comfort from heel to toe.

## Vectra Pro Options:

Available in tall or short designs with your choice of air or non-air models, and either aluminum or plastic struts for tailored support.

AIR ALUMINUM STRUT



NON-AIR ALUMINUM STRUT



AIR PLASTIC STRUT



NON-AIR PLASTIC STRUT



**1** SIDE-MOUNTED AIR PUMP AND RELEASE VALVE allows effortless adjustments on pneumatic models

**2** MEDIAL AND LATERAL AIR BLADDERS with foam lining ensure comfortable, enhanced stabilization

**3** ALUMINUM AND PLASTIC STRUTS offer customizable support

**4** HARD TOE COVER extends protection and blocks debris



**5** REMOVABLE HEEL CLIP offers customizable fit and support while accommodating compatible Achilles wedges, sold separately



**6** WIDE FOOTBED comfortably accommodates swelling and bandages

**7** ROCKER BOTTOM promotes natural gait for improved mobility







Breg, Inc, 2382 Faraday Ave, Ste # 300, Carlsbad, CA 92008, USA  
P: 800.321.0607, F: 760.795.5295, E: orderprocess@BREG.com, W: BREG.com



## ORDERING INFORMATION

BL1000-0X0	Short, Non-Air, Plastic Strut
BL2000-0X0	Tall, Non-Air, Plastic Strut
BL3000-0X0	Short, Air, Plastic Strut
BL4000-0X0	Tall, Air, Plastic Strut
BL5000-0X0	Short, Non-Air, Aluminum Strut
BL6000-0X0	Tall, Non-Air, Aluminum Strut
BL7000-0X0	Short, Air, Aluminum Strut
BL8000-0X0	Tall, Air, Aluminum Strut

## ACCESSORIES

BLO001-0X0	Short Boot Replacement Kit
BLO002-0X0	Tall Boot Replacement Kit
BLO005-0X0	Extension Plate (not available in XS)

## SIZE CHART

PART X	BOOT SIZE	MEN'S	WOMEN'S	EU
1	XS	< 4	< 5.5	35 - 36
2	S	4.5 - 7	6 - 8	37 - 40
3	M	7.5 - 10	8.5 - 11.5	40 - 43
4	L	10.5 - 12.5	12 - 13.5	43 - 46
5	XL	13 +	14 +	47 +

Learn More:



PDAC Available

©2025 BREG, Inc. All rights reserved. 03/25, Vectra is registered trademark of BREG, Inc.  
P/N: 10.000017 Rev B



MDSS GmbH  
Schiffgraben 41  
D-30175 Hannover  
Germany