



Breg, Inc, 2382 Faraday Ave, Ste # 300,
 Carlsbad, CA 92008, USA
 P: 800.321.0607, F: 760.795.5295,
 E: orderprocess@BREG.com, W: BREG.com



Skyrocket to Recovery

Designed for little feet, built for big adventures—helping kids recover with comfort, style, and durability.

VECTRA WEE WALKER BOOTS

BL0100-020	Small Vectra Wee Walker Boot
BL0100-030	Medium Vectra Wee Walker Boot
BL0100-040	Large Vectra Wee Walker Boot
BL0200-020	AIR Small Vectra Wee Walker Boot
BL0200-030	AIR Medium Vectra Wee Walker Boot
BL0200-040	AIR Large Vectra Wee Walker Boot

REPLACEMENT KITS

BL0006-020	Small Replacement Kit
BL0006-030	Medium Replacement Kit
BL0006-040	Large Replacement Kit
BL0007-020	AIR Small Replacement Kit
BL0007-030	AIR Medium Replacement Kit
BL0007-040	AIR Large Replacement Kit

PATCH ACCESSORIES

BL0000-100	Pink Patches - Set of 3
BL0000-200	Gray Patches - Set of 3

SIZE CHART

BOOT SIZE	FOOT LENGTH
S	Up to 6 inches (15.2 cm)
M	6 - 7 inches (15.2 - 17.8 cm)
L	7 - 8 inches (17.8 - 20.3 cm)



MDSS GmbH
 Schiffgraben 41
 D-30175 Hannover
 Germany



Vectra® Wee Walker

Introducing **Vectra Wee Walker**

The Vectra Wee Walker is a pediatric walking boot designed for both comfort and durability. Built to support active children, it combines sturdy engineering to withstand the wear and tear of everyday play. Its adjustable fit ensures a secure, comfortable experience, to help with recovery while staying active. Available with fun, space-themed patches, it's the perfect companion to keep them moving forward in their recovery.

1 SIDE-MOUNTED AIR PUMP AND RELEASE VALVE allows effortless adjustments on pneumatic models

2 LIGHTWEIGHT DESIGN ensures ease of movement while withstanding the wear and tear of everyday play

3 HOOK & LOOP STRAPS for a secure, customizable fit

4 DURABLE NON-SLIP TREAD for enhanced traction and stability on various surfaces

5 ROCKER BOTTOM promotes natural gait for improved mobility



Vectra Wee Walker Options:

Available in both air and non-air models with size specific designs: the small with a spaceship design, the medium with an astronaut design, and the large with a rocketship design.



SMALL AIR AND NON-AIR




MEDIUM AIR AND NON-AIR



LARGE AIR AND NON-AIR

Vectra Wee Walker Patches:

Add a little extra fun with optional patches! Choose from the pink pack of 3 or the gray pack of 3, each featuring fun space-themed designs to give each walker boot a personalized touch.



PINK PATCHES



GRAY PATCHES